

Product Datasheet

DAPK3/ZIPK Antibody NBP1-03077

Unit Size: 0.05 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-03077

Updated 2/24/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-03077



NBP1-03077**DAPK3/ZIPK Antibody**


Product Information	
Unit Size	0.05 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS

Product Description	
Host	Rabbit
Gene ID	1613
Gene Symbol	DAPK3
Species	Human, Mouse, Rat
Immunogen	Synthetic peptide from human DAPK3/ZIPK. The immunogen is located within aa 270-320 of DAPK3/ZIPK.


Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunocytochemistry
Recommended Dilutions	Western Blot 1:500-1:1000, ELISA, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:10-1:500, Immunohistochemistry-Paraffin 2.5 ug/ml, Immunocytochemistry

Images


Immunocytochemistry/Immunofluorescence: DAPK3/ZIPK Antibody [NBP1-03077] - Immunofluorescence of DAPK3/ZIPK in Jurkat cells with DAPK3/ZIPK Antibody at 10 ug/ml.




Immunohistochemistry-Paraffin: DAPK3/ZIPK Antibody [NBP1-03077] - Staining of human brain, cortex. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentration 5 ug/ml.




Immunocytochemistry: DAPK3/ZIPK Antibody [NBP1-03077] - Immunocytochemistry of DAPK3/ZIPK in Jurkat cells with ZIPK antibody at 10 ug/mL.



Immunocytochemistry/Immunofluorescence: DAPK3/ZIPK Antibody [NBP1-03077] - Immunofluorescence of DAPK3/ZIPK in Jurkat cells with DAPK3/ZIPK antibody at 5 ug/ml. Green: . Blue: DAPI staining.



Immunohistochemistry-Paraffin: DAPK3/ZIPK Antibody [NBP1-03077] - Staining of human small intestine, Peyer's patch. Immunohistochemistry of formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentration 5 ug/ml.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-03077

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

