

# Product Datasheet

## Geminin Antibody (2H7) - BSA Free NBP1-04293

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-04293](http://www.novusbio.com/NBP1-04293)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-04293](http://www.novusbio.com/reviews/destination/NBP1-04293)



**NBP1-04293**

Geminin Antibody (2H7) - BSA Free


| Product Information |                                                                                        |
|---------------------|----------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                 |
| Concentration       | 1.0 mg/ml                                                                              |
| Storage             | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality           | Monoclonal                                                                             |
| Clone               | 2H7                                                                                    |
| Preservative        | 0.02% Sodium Azide                                                                     |
| Isotype             | IgG1 Kappa                                                                             |
| Purity              | Protein G purified                                                                     |
| Buffer              | PBS (pH 7.4), 10% Glycerol                                                             |

| Product Description |                                                            |
|---------------------|------------------------------------------------------------|
| Host                | Mouse                                                      |
| Gene ID             | 51053                                                      |
| Gene Symbol         | GMNN                                                       |
| Species             | Human                                                      |
| Immunogen           | Recombinant human Geminin (1-209 aa) purified from E. coli |

| Product Application Details |                                               |
|-----------------------------|-----------------------------------------------|
| Applications                | Western Blot, ELISA                           |
| Recommended Dilutions       | Western Blot 1:500-1:1000, ELISA 1:100-1:2000 |

**Images**

Western Blot: Geminin Antibody (2H7) [NBP1-04293] - The 293T cell lysates (50ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human Geminin (1:500). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-04293**

---

|                  |                                                         |
|------------------|---------------------------------------------------------|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B          | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-43319-0.5mg | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)           |
| NBC1-18423       | Recombinant Human Geminin His Protein                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-04293](http://www.novusbio.com/reviews/submit/NBP1-04293)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

