

# Product Datasheet

## UbcH8/Ube2L6 Antibody (k1H3) - BSA Free NBP1-04351

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-04351](http://www.novusbio.com/NBP1-04351)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-04351](http://www.novusbio.com/reviews/destination/NBP1-04351)



**NBP1-04351**

UbcH8/Ube2L6 Antibody (k1H3) - BSA Free

**Product Information**

<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	k1H3
<b>Preservative</b>	0.02% Sodium Azide
<b>Isotype</b>	IgG2b Kappa
<b>Purity</b>	Protein G purified
<b>Buffer</b>	PBS (pH 7.4), 10% Glycerol

**Product Description**


<b>Host</b>	Mouse
<b>Gene ID</b>	9246
<b>Gene Symbol</b>	UBE2L6
<b>Species</b>	Human
<b>Immunogen</b>	Recombinant human UbcH8/Ube2L6 (1-152aa) purified from E. coli

**Product Application Details**


<b>Applications</b>	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
<b>Recommended Dilutions</b>	Western Blot 1:2000, ELISA 1:100-1:2000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin 1:10-1:500

**Images**


Western Blot: UbcH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - The Cell lysate (20ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human UBE2L6 antibody (1:2000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: Jurkat cell lysate




**Immunocytochemistry/Immunofluorescence:** UbchH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - ICC/IF analysis of UBE2L6 in HeLa cells line, stained with DAPI (Blue) for nucleus staining and monoclonal anti-human UBE2L6 antibody (1:100) with goat anti-mouse IgG-Alexa fluor 488 conjugate (Green).




**Immunohistochemistry-Paraffin:** UbchH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - Ube2L6 Antibody (k1H3) [NBP1-04351]-Paraffin embedded sections of human breast lobule tissue were incubated with anti-human UBE2L6 for 2 hours at room temperature. Antigen retrieval was performed in 0.1M sodium citrate buffer and detected using Diaminobenzidine (DAB)



**Western Blot:** UbchH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - Ube2L6 Antibody (k1H3) [NBP1-04351]-Cell lysates of MCF7 (20ug) were resolved by SDS-PAGE, transferred to NC membrane and probed with anti-human UBE2L6. Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.



**Western Blot:** UbchH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - The Cell lysates (5ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human UBE2L6 antibody (1:2000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: 293T cell lysate Lane 2.: UBE2L6 Transfected 293T cell lysate





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-04351**

---

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
NBP1-86885PEP	Ubch8/Ube2L6 Recombinant Protein Antigen

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-04351](http://www.novusbio.com/reviews/submit/NBP1-04351)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

