

# Product Datasheet

## HPi4 Antibody (HIC1-5F10) [DyLight 488] NBP1-18948G

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-18948G](http://www.novusbio.com/NBP1-18948G)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-18948G](http://www.novusbio.com/reviews/destination/NBP1-18948G)



**NBP1-18948G**

HPi4 Antibody (HIC1-5F10) [DyLight 488]

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	HIC1-5F10
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 488
Purity	Tissue culture supernatant
Buffer	50mM Sodium Borate
Product Description	
Description	This product is provided as conjugated tissue culture supernatant. Due to the crude preparation of the antibody material, excess dye may remain in the final product. We highly recommend taking necessary precautions, such as adding additional wash steps and running proper controls.
Host	Mouse
Species	Human, Mouse
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 26515645).
Marker	Endocrine pan-islet cell Marker
Specificity/Sensitivity	Specific for multiple endocrine cell types.
Immunogen	Human pancreatic enriched islet cells containing low levels of exocrine and ductal cells.
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Frozen
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-18948G**

---

NBP1-97005G	Mouse IgG1 Isotype Control (MG1) [DyLight 488]
NBP1-18948R	HPi4 Antibody (HIC1-5F10) [DyLight 550]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-18948G](http://www.novusbio.com/reviews/submit/NBP1-18948G)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

