

Product Datasheet

HPa2 Antibody (DHIC2-2B4) [DyLight 350] NBP1-18950UV

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-18950UV

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-18950UV



NBP1-18950UV

HPa2 Antibody (DHIC2-2B4) [DyLight 350]

| Product Information | |
|---------------------|--|
| Unit Size | 0.1 ml |
| Concentration | This product is unpurified. The exact concentration of antibody is not quantifiable. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | DHIC2-2B4 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | DyLight 350 |
| Purity | Tissue culture supernatant |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|-------------------------|---|
| Description | This product is provided as conjugated tissue culture supernatant. Due to the crude preparation of the antibody material, excess dye may remain in the final product. We highly recommend taking necessary precautions, such as adding additional wash steps and running proper controls. |
| Host | Mouse |
| Species | Human |
| Marker | Alpha Cell Marker |
| Specificity/Sensitivity | Specific for alpha endocrine cell types. |
| Immunogen | Human pancreatic enriched alpha cells |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen |
| Recommended Dilutions | Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Frozen |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-18950UV

| | |
|--------------|--|
| NBP1-97005UV | Mouse IgG1 Isotype Control (MG1) [DyLight 350] |
| NBP1-18950B | HPa2 Antibody (DHIC2-2B4) [Biotin] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-18950UV

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

