

# Product Datasheet

## HPd2 Antibody (DHIC2-5H10) [DyLight 350] NBP1-18955UV

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-18955UV](http://www.novusbio.com/NBP1-18955UV)

Updated 9/10/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-18955UV](http://www.novusbio.com/reviews/destination/NBP1-18955UV)



**NBP1-18955UV**

HPd2 Antibody (DHIC2-5H10) [DyLight 350]

| Product Information         |                                                                                                                                                                                                                                                                                           |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                    |
| Concentration               | This product is unpurified. The exact concentration of antibody is not quantifiable.                                                                                                                                                                                                      |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                  |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                |
| Clone                       | DHIC2-5H10                                                                                                                                                                                                                                                                                |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                        |
| Isotype                     | IgM                                                                                                                                                                                                                                                                                       |
| Conjugate                   | DyLight 350                                                                                                                                                                                                                                                                               |
| Purity                      | Tissue culture supernatant                                                                                                                                                                                                                                                                |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                                                                                                                                           |
| Description                 | This product is provided as conjugated tissue culture supernatant. Due to the crude preparation of the antibody material, excess dye may remain in the final product. We highly recommend taking necessary precautions, such as adding additional wash steps and running proper controls. |
| Host                        | Mouse                                                                                                                                                                                                                                                                                     |
| Species                     | Human                                                                                                                                                                                                                                                                                     |
| Reactivity Notes            | Human.                                                                                                                                                                                                                                                                                    |
| Marker                      | Pancreatic, liver duct cell Marker                                                                                                                                                                                                                                                        |
| Specificity/Sensitivity     | Specific for duct cell types.                                                                                                                                                                                                                                                             |
| Immunogen                   | Human pancreatic enriched duct cells                                                                                                                                                                                                                                                      |
| Notes                       | Dylight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                         |
| Product Application Details |                                                                                                                                                                                                                                                                                           |
| Applications                | Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen                                                                                                                                                                                 |
| Recommended Dilutions       | Flow Cytometry, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Frozen                                                                                                                                                                                 |
| Application Notes           | This HPd2 (DHIC2-5H10) antibody is useful for Immunohistochemistry on acetone fixed frozen sections, Immunocytochemistry/Immunofluorescence and Flow cytometry.                                                                                                                           |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-18955UV](http://www.novusbio.com/reviews/submit/NBP1-18955UV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

