

Product Datasheet

CD117/c-kit Antibody (kit2c75) [FITC] NBP1-28328

Unit Size: 0.25 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-28328

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-28328



NBP1-28328

CD117/c-kit Antibody (kit2c75) [FITC]

Product Information

| | |
|----------------------|---|
| Unit Size | 0.25 mg |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Monoclonal |
| Clone | kit2c75 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a Kappa |
| Conjugate | FITC |
| Purity | Protein A or G purified |
| Buffer | PBS |

Product Description


| | |
|--------------------------------|---|
| Host | Mouse |
| Gene ID | 3815 |
| Gene Symbol | KIT |
| Species | Chicken |
| Marker | Hematopoietic Stem Cell Marker |
| Specificity/Sensitivity | Chicken c-Kit |
| Immunogen | Sp/chkit4A4 cells (Sp2/0 cells transfected with chicken c-Kit cDNA) |

Product Application Details

| | |
|------------------------------|-------------------------------------|
| Applications | Flow Cytometry, Knockdown Validated |
| Recommended Dilutions | Flow Cytometry, Knockdown Validated |

Images

CD117/c-kit Antibody (kit2c75) [FITC] [NBP1-28328] - Vial of FITC conjugated antibody. FITC is optimally excited at 498 nm by the Blue laser (488 nm) and has an emission maximum of 519 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-28328

| | |
|---------------|---|
| NBP1-96981F | Mouse IgG2a Kappa Isotype Control (M2AK) [FITC] |
| NBP1-28329 | CD117/c-kit Antibody (kit2c75) [Biotin] |
| NBP2-52974PEP | CD117/c-kit Recombinant Protein Antigen |
| 210-TA-005 | TNF-alpha [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-28328

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

