

Product Datasheet

Respiratory Syncytial Virus Antibody - BSA Free NBP1-30413

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-30413

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-30413



NBP1-30413**Respiratory Syncytial Virus Antibody - BSA Free**


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Protein G purified
Buffer	PBS

Product Description	
Host	Rabbit
Species	Virus
Reactivity Notes	Virus
Specificity/Sensitivity	Whole RSV virions (subgroup A-Long strain) were used as an immunogen. Detects G, F, N and M2 proteins of RSV.
Immunogen	Whole RSV virions (subgroup A-Long strain) were used as an immunogen.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:500 -1:1000
Application Notes	This antibody is useful in Western blot.

Images

Western Blot: Respiratory Syncytial Virus Antibody [NBP1-30413] - For WB analysis of Respiratory Syncytial Virus (RSV), 3 ug of RSV virions (subgroup A-Long strain) were loaded on SDS-PAGE, transferred to a PVDF membrane and detected with NBP1-30413 at a dilution of 1:1000. Goat anti-rabbit secondary antibody was used at 1:10,000 dilution. The antibody detects the following RSV proteins: Glycoprotein (G) ~90 kDa, Fusion (F) protein) ~55 kDa, and nucleocapsid (N) protein) ~46 kDa. Overexposure of the blot can result in detection of M2 protein (~22 kDa-not shown).

**Publications**

Garofalo R, Haeberle H. Epithelial regulation of innate immunity to respiratory syncytial virus. Am J Respir Cell Mol Biol. 2000-01-01 [PMID: 11062135]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-30413

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-30413

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

