

# Product Datasheet

## UbcH7/UBE2L3 Antibody - BSA Free NBP1-30995

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-30995](http://www.novusbio.com/NBP1-30995)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-30995](http://www.novusbio.com/reviews/destination/NBP1-30995)



**NBP1-30995**

UbchH7/UBE2L3 Antibody - BSA Free

**Product Information**

|                                |                                                                                                                  |
|--------------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>               | 100 ul                                                                                                           |
| <b>Concentration</b>           | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| <b>Clonality</b>               | Polyclonal                                                                                                       |
| <b>Preservative</b>            | 0.01% Thimerosal                                                                                                 |
| <b>Isotype</b>                 | IgG                                                                                                              |
| <b>Purity</b>                  | Antigen Affinity-purified                                                                                        |
| <b>Buffer</b>                  | PBS (pH 7), 40% Glycerol                                                                                         |
| <b>Target Molecular Weight</b> | 18 kDa                                                                                                           |

**Product Description**


|                         |                                                                                                                                                                                                                   |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Novus Biologicals Rabbit UbchH7/UBE2L3 Antibody - BSA Free (NBP1-30995) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>             | Rabbit                                                                                                                                                                                                            |
| <b>Gene ID</b>          | 7332                                                                                                                                                                                                              |
| <b>Gene Symbol</b>      | UBE2L3                                                                                                                                                                                                            |
| <b>Species</b>          | Human, Mouse, Zebrafish                                                                                                                                                                                           |
| <b>Reactivity Notes</b> | Immunogen displays the following percentage of sequence identity for non-tested species: <i>Xenopus laevis</i> (100%) and <i>Xenopus Tropicalis</i> (100%).                                                       |
| <b>Immunogen</b>        | Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human UbchH7/UBE2L3. The exact sequence is proprietary.                                                      |

**Product Application Details**


|                              |                                                                                                                                                                                                 |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation                                                                  |
| <b>Recommended Dilutions</b> | Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunoprecipitation 1:100-1:500, Immunohistochemistry-Paraffin 1:100-1:1000 |

**Images**


Western Blot: UbchH7/UBE2L3 Antibody [NBP1-30995] - Mouse tissue extract (50 ug) was separated by 15% SDS-PAGE, and the membrane was blotted with UbchH7/UBE2L3 antibody (GTX104717).




Immunocytochemistry/Immunofluorescence: Ubch7/UBE2L3 Antibody [NBP1-30995] - Analysis of HeLa, using NBP1-30995 at 1:200 dilution.




Immunohistochemistry-Paraffin: Ubch7/UBE2L3 Antibody [NBP1-30995] - A549 xenograft, using UBE2L3 antibody at 1:100 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Western Blot: Ubch7/UBE2L3 Antibody [NBP1-30995] - A. 30 ug Neuro2A whole cell lysate/extract, B. 30 ug GL261 whole cell lysate/extract, C. 30 ug C8D30 whole cell lysate/extract 15 % SDS-PAGE gel, antibody dilution 1:1000.




Western Blot: Ubch7/UBE2L3 Antibody [NBP1-30995] - Sample (30 ug of whole cell lysate) A: HeLa B: Hep G2 12% SDS PAGE; antibody diluted at 1:1000.



Western Blot: UbchH7/UBE2L3 Antibody [NBP1-30995] - Whole zebrafish extract (30 ug) was separated by 15% SDS-PAGE, and the membrane was blotted with UbchH7/UBE2L3 antibody (GTx104717).



Immunoprecipitation: UbchH7/UBE2L3 Antibody [NBP1-30995] - Immunoprecipitation of UbchH7/UBE2L3 protein from HeLa whole cell extracts using 5 ug of UbchH7/UBE2L3 antibody (NBP1-30995). Western blot analysis was performed using UbchH7/UBE2L3 antibody (NBP1-30995). EasyBlot anti-Rabbit IgG was used as a secondary reagent.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-30995**

---

|                 |                                                      |
|-----------------|------------------------------------------------------|
| HAF008          | Goat anti-Rabbit IgG Secondary Antibody [HRP]        |
| NB7160          | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]  |
| NBP2-24891      | Rabbit IgG Isotype Control                           |
| NBP2-35110-10ug | Recombinant Human UbchH7/UBE2L3 His (N-Term) Protein |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-30995](http://www.novusbio.com/reviews/submit/NBP1-30995)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

