

Product Datasheet

IFIT2 Antibody - BSA Free NBP1-31164

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 2 Publications: 6

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-31164

Updated 2/21/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-31164



NBP1-31164


IFIT2 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris-Glycine (pH7), 10% Glycerol
Target Molecular Weight	55 kDa
Product Description	
Host	Rabbit
Gene ID	3433
Gene Symbol	IFIT2
Species	Human
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the N-terminus region of human IFIT2. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500-1:3000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000




Images


Western Blot: IFIT2 Antibody [NBP1-31164] - p115 and Sec16A are required for optimal IRF-3 activation in response to activation of cytosolic RNA and DNA sensors. HeLa cells were infected with lentiviral vectors encoding shRNA targeting p115. Image collected and cropped by CiteAb from the following publication (<https://doi.org/10.1371/journal.ppat.1002747>) licensed under a CC-BY license.




Immunocytochemistry/Immunofluorescence: IFIT2 Antibody [NBP1-31164] - HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: IFIT2 protein stained by IFIT2 antibody [N1], N-term diluted at 1:200. Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody diluted at 1:10000. Blue: Hoechst 33342 staining.



Western Blot: IFIT2 Antibody [NBP1-31164] - Sample(30 ug of whole cell lysate)A:A431 B:Hep G2 7.5% SDS PAGE, antibody diluted at 1:1000.



Western Blot: IFIT2 Antibody [NBP1-31164] - p115 & Sec16A are required for optimal IRF-3 activation in response to activation of cytosolic RNA & DNA sensors. HeLa cells were infected with lentiviral vectors encoding shRNA targeting p115 (A) or Sec16A (B) & non-targeting (NT) control shRNA for 24 h followed by puromycin selection (1.5 μ g/ml) for 5 days. Cells were left untreated or stimulated with poly I:C (1 μ g/ml), poly dA:dT (1 μ g/ml) or SeV (200 HAU/ml) for 16 h. Whole-cell lysates were prepared & subjected to immunoblot analysis with indicated antibodies. One of two independent experiments with similar results is shown. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/22792062/>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Rex EA, Seo D, Chappidi S, Pinkham C et Al. FEAR antiviral response pathway is independent of interferons and countered by poxvirus proteins Nat Microbiol 2024-03-28 [PMID: 38538832]

Khowal S, Monga S, Naqvi S et al. Molecular winnowing, expressional analyses and interactome scrutiny of cellular proteomes of oral squamous cell carcinoma Advances in Cancer Biology - Metastasis 2021-10-01

Khan KA TRK-Fused Gene (TFG), a protein involved in protein secretion pathways, is an essential component of the antiviral innate immune response PloS Pathog 2021-01-07 [PMID: 33411856] (WB, WB, Human)

Khan KA, Do F, Marineau A et al. Fine-Tuning of the RIG-I-Like Receptor/Interferon Regulatory Factor 3-Dependent Antiviral Innate Immune Response by the Glycogen Synthase Kinase 3/beta-Catenin Pathway. Mol Cell Biol 2015-09-01 [PMID: 26100021]

Baril M, Es-Saad S, Chatel-Chaix L et al. Genome-wide RNAi Screen Reveals a New Role of a WNT/CTNNB1 Signaling Pathway as Negative Regulator of Virus-induced Innate Immune Responses. PLoS Pathog 2013-06-01 [PMID: 23785285] (WB, Human)

van Zuylen WJ, Doyon P, Clement JF et al. Proteomic profiling of the TRAF3 interactome network reveals a new role for the ER-to-Golgi transport compartments in innate immunity. PLoS Pathog 2012-01-01 [PMID: 22792062] (WB, Human)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-31164

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00003433-P01-10ug	Recombinant Human IFIT2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-31164

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

