

Product Datasheet

NQO-1 Antibody - BSA Free

NBP1-31355

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-31355

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-31355



NBP1-31355

NQO-1 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol
Target Molecular Weight	31 kDa

Product Description


Description	Novus Biologicals Rabbit NQO-1 Antibody - BSA Free (NBP1-31355) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-NQO-1 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1728
Gene Symbol	NQO1
Species	Human, Mouse, Rat, Monkey
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human NQO-1. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Knockdown Validated

Images


Immunohistochemistry: NQO-1 Antibody [NBP1-31355] - (b) Immunohistochemical analyses of RKO/pNQO1, RKO/shCont, and RKO/pshNQO1 xenograft tumours. The sections were stained for NQO1. Scale bar, 50 μ m. (c) Quantification of NQO1 in RKO/pNQO1, RKO/shCont and RKO/pshNQO1 xenograft tumours (n=3 each group). n=5 in each tumour. Two-tail t-test. *P<0.05. All error bars represent the mean \pm s.e.m. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/ncomms13593>) licensed under a CC-BY license.




Western Blot: NQO-1 Antibody [NBP1-31355] - Various whole cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with NQO1 antibody [C2C3], C-term diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Immunohistochemistry: NQO-1 Antibody [NBP1-31355] - IHC detection of high-level expression of NQO1 (n=11) compared with the low-level expression of NQO1 (n=9). Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/ncomms13593>) licensed under a CC-BY license.




Immunocytochemistry/Immunofluorescence: NQO-1 Antibody [NBP1-31355] - Paraformaldehyde-fixed HeLa, using antibody at 1:200 dilution.




Western Blot: NQO-1 Antibody [NBP1-31355] - A. 50 ug mouse heart lysate/extract 12% SDS-PAGE NQO1 antibody [C2C3], C-term dilution: 1:500 The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Western Blot: NQO-1 Antibody [NBP1-31355] - A. 50 ug rat heart lysate/extract 12% SDS-PAGE NQO1 antibody [C2C3], C-term (GTX100235) dilution: 1:500 The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Immunohistochemistry-Paraffin: NQO-1 Antibody [NBP1-31355] - Human rectal carcinoma. NQO1 antibody [C2C3], C-term diluted at 1:500. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Immunohistochemistry-Paraffin: NQO-1 Antibody [NBP1-31355] - Mouse kidney. NQO1 stained by NQO1 antibody [C2C3], C-term diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min.



Immunohistochemistry-Paraffin: NQO-1 Antibody [NBP1-31355] - Mouse kidney. NQO1 stained by NQO1 antibody [C2C3], C-term diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min.



Immunohistochemistry: NQO-1 Antibody [NBP1-31355] - NQO1 promotes in vivo tumour growth. (a) RKO/pNQO1, RKO/shCont & RKO/pshNQO1 cells injected subcutaneously into right flank of athymic, 7-week-old female BALB/C nude mice, & tumour growth assessed. Tumour volume (TV) calculated by using following formula: $TV = \text{length} \times (\text{width})^2 \times 0.5$. Each group contained 12 animals. (** $P < 0.01$ w/ unpaired t-test). (b) Immunohistochemical analyses of RKO/pNQO1, RKO/shCont, & RKO/pshNQO1 xenograft tumours. The sections stained for NQO1, HIF-1 α , proliferation (Ki67) & apoptosis (CC3) using 3,3'-DAB. Scale bar, 50 μm . (c) Quantification of NQO1, proliferative marker Ki67 & apoptotic marker CC3 in RKO/pNQO1, RKO/shCont & RKO/pshNQO1 xenograft tumours ($n=3$ each group). $n=5$ in each tumour. Two-tail t-test. ** $P < 0.01$, * $P < 0.05$. #, not significant. All error bars represent mean \pm s.e.m. (d) RKO/pshCont1/pshCont2, RKO/pNQO1/pshCont2 & RKO/pNQO1/pshHIF-1 α injected subcutaneously into right flank of athymic, 7-week-old female BALB/C nude mice, & tumour growth assessed. Tumour volume (TV) calculated by using following formula: $TV = \text{length} \times (\text{width})^2 \times 0.5$. Each group contained 10 animals (* $P < 0.05$ w/ unpaired t-test). (e) Immunohistochemical analyses of RKO/pshCont1/pshCont2, RKO/pNQO1/pshCont2 & RKO/pNQO1/pshHIF-1 α xenograft tumours for HIF-1 α , proliferation (Ki67), apoptosis (CC3) & vasculature (CD34) using 3,3'-DAB. Arrowheads denote blood vessels. Scale bar, 50 μm . (f) RKO/pNQO1/pshCont or RKO/pNQO1/pshHIF-1 α xenograft tumours ($n=3-5$ each group) quantified for proliferation (Ki67), apoptosis (CC3) & vascularization (CD34). $n=5$ in each tumour. Two tailed t-test. ** $P < 0.01$, * $P < 0.05$. #, not significant. All error bars represent mean \pm s.e.m. (g) Schematic model showing how NQO1 stabilizes HIF-1 α in cancer cells. Image collected & cropped by CiteAb from following publication (<https://pubmed.ncbi.nlm.nih.gov/27966538>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Oh ET, Kim HG, Kim CH et al. NQO1 regulates cell cycle progression at the G2/M phase Theranostics 2023-01-10 [PMID: 36793872] (Immunohistochemistry-Paraffin, Human)

Oh ET, Kim JW, Kim JM et al. NQO1 inhibits proteasome-mediated degradation of HIF-1 α . Nat Commun. 2016-12-14 [PMID: 27966538] (ICC/IF, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-31355

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-85223PEP	NQO-1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-31355

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

