

Product Datasheet

F-box protein 43 Antibody NBP1-31473

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-31473

Updated 11/30/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-31473



NBP1-31473

F-box protein 43 Antibody

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7), 10% Glycerol
Target Molecular Weight	78 kDa

Product Description


Host	Rabbit
Gene ID	286151
Gene Symbol	FBXO43
Species	Human, Mouse
Reactivity Notes	Xenopus laevis (100%).
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human FBXO43. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Western Blot: F-box protein 43 Antibody [NBP1-31473] - Sample (30 ug of whole cell lysate) A: 293T 7. 5% SDS PAGE; antibody diluted at 1:1000.




Immunocytochemistry/Immunofluorescence: F-box protein 43 Antibody [NBP1-31473] - Paraformaldehyde-fixed A431, using antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: F-box protein 43 Antibody [NBP1-31473] - Paraffin-embedded TOV-21G xenograft. FBXO43 antibody [C3], C-term diluted at 1:250.



Western Blot: F-box protein 43 Antibody [NBP1-31473] - Sample (30 ug of whole cell lysate) A: NIH-3T3 B: JC 7. 5% SDS PAGE; antibody diluted at 1:1000.



Immunohistochemistry-Paraffin: F-box protein 43 Antibody [NBP1-31473] - Paraffin-embedded Saos-2 xenograft, using antibody at 1:100 dilution.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-31473

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

