

Product Datasheet

GNAQ Antibody - BSA Free

NBP1-31985

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-31985

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-31985



NBP1-31985

GNAQ Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol
Target Molecular Weight	42 kDa

Product Description


Description	Novus Biologicals Rabbit GNAQ Antibody - BSA Free (NBP1-31985) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	2776
Gene Symbol	GNAQ
Species	Human, Mouse, Rat
Reactivity Notes	Xenopus laevis (96%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human GNAQ. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Western Blot: GNAQ Antibody [NBP1-31985] - A. 30 ug 293T whole cell extract B. 30 ug whole cell extract of 3xFlag-human GNAQ-transfected 293T cells 10 % SDS-PAGE GNAQ antibody dilution: 1:1000




Immunohistochemistry-Paraffin: GNAQ Antibody [NBP1-31985] - Paraffin-embedded mouse ovary. GNAQ antibody dilution: 1:500.




Western Blot: GNAQ Antibody [NBP1-31985] - Whole cell extract (30 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with GNAQ antibody diluted at 1:500.




Immunohistochemistry-Paraffin: GNAQ Antibody [NBP1-31985] - Paraffin-embedded placenta. GNAQ antibody dilution: 1:100.



Immunohistochemistry-Paraffin: GNAQ Antibody [NBP1-31985] - Paraffin-embedded rat intestine. GNAQ antibody dilution: 1:500.



Immunohistochemistry-Paraffin: GNAQ Antibody [NBP1-31985] -
Paraffin-embedded mouse prostate. GNAQ antibody dilution: 1:500.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-31985

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-51808-0.1mg	Recombinant Human GNAQ His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-31985

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

