

Product Datasheet

MAD1L1/MAD1 Antibody NBP1-32775

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-32775

Updated 8/8/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-32775



NBP1-32775**MAD1L1/MAD1 Antibody**


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7), 20% Glycerol
Target Molecular Weight	83 kDa

Product Description	
Host	Rabbit
Gene ID	8379
Gene Symbol	MAD1L1
Species	Human, Mouse, Rat
Reactivity Notes	Potorous tridactylus/Potoroo (100%), Bovine (86%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human MAD1. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Western Blot 1:500-1:20000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/Immunofluorescence 1:100-1:1000, Immunoprecipitation 1:100-1:500, Immunohistochemistry-Paraffin 1:100-1:1000, Chromatin Immunoprecipitation (ChIP) 1:10-1:500
Application Notes	ChIP assay- Assay dependent

Images


Western Blot: MAD1L1/MAD1 Antibody [NBP1-32775] - Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with MAD1 antibody diluted by 1:5000.




Immunocytochemistry/Immunofluorescence: MAD1L1/MAD1 Antibody [NBP1-32775] - HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: MAD1 protein stained by MAD1 antibody diluted at 1:500. Red: HEC1, a kinetochore marker, stained by HEC1 antibody [9G3.23] diluted at 1:500. Blue: Hoechst 33342 staining.




Immunohistochemistry-Paraffin: MAD1L1/MAD1 Antibody [NBP1-32775] - Paraffin-embedded SW480 Xenograft, using antibody at 1:100 dilution.




Western Blot: MAD1L1/MAD1 Antibody [NBP1-32775] - Sample (30 ug of whole cell lysate) A: Molt-4 7.5% SDS PAGE diluted at 1:10000




Western Blot: MAD1L1/MAD1 Antibody [NBP1-32775] - A. 30 ug PC-12 whole cell lysate/extract B. 30 ug Rat 2 whole cell lysate/extract 7.5% SDS-PAGE MAD1 antibody dilution: 1:3000 The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Western Blot: MAD1L1/MAD1 Antibody [NBP1-32775] - A. 30 ug Neuro2A whole cell lysate/extract B. 30 ug GL261 whole cell lysate/extract C. 30 ug C8D30 whole cell lysate/extract D. 30 ug NIH-3T3 whole cell lysate/extract E. 30 ug BCL-1 whole cell lysate/extract F. 30 ug Raw264.7 whole cell lysate/extract 7.5% SDS-PAGE MAD1 antibody dilution: 1:3000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Immunoprecipitation: MAD1L1/MAD1 Antibody [NBP1-32775] - Immunoprecipitation of MAD1 protein from HeLa whole cell extracts using 5 ug of MAD1 antibody. Western blot analysis was performed using MAD1 antibody. EasyBlot anti-Rabbit IgG was used as a secondary reagent.



Publications

Chen Y, Yeh PC, Huang JC et al. The spindle checkpoint protein MAD1 regulates the expression of E-cadherin and prevents cell migration. *Oncol Rep* 2012-02-01 [PMID: 22024768] (Chemotaxis, WB, ICC/IF, Human)

Yeh PC, Yeh CC, Chen YC, Juang YL. RED, a spindle pole-associated protein, is required for kinetochore localization of MAD1, mitotic progression, and activation of the spindle assembly checkpoint. *J Biol Chem*. 2012-02-18 [PMID: 22351768] (ICC/IF, Chemotaxis, WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-32775

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00008379-Q01-10ug	Recombinant Human MAD1L1/MAD1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-32775

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

