

# Product Datasheet

## SOX13 Antibody - BSA Free

### NBP1-32847

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-32847](http://www.novusbio.com/NBP1-32847)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-32847](http://www.novusbio.com/reviews/destination/NBP1-32847)



**NBP1-32847**


SOX13 Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                      |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                                                                               |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                                                     |
| Storage                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                                                                                                         |
| Clonality                   | Polyclonal                                                                                                                                                                                           |
| Preservative                | 0.01% Thimerosal                                                                                                                                                                                     |
| Isotype                     | IgG                                                                                                                                                                                                  |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                            |
| Buffer                      | 0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol                                                                                                                                                         |
| Target Molecular Weight     | 69 kDa                                                                                                                                                                                               |
| Product Description         |                                                                                                                                                                                                      |
| Description                 | Novus Biologicals Rabbit SOX13 Antibody - BSA Free (NBP1-32847) is a polyclonal antibody validated for use in WB, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                        | Rabbit                                                                                                                                                                                               |
| Gene ID                     | 9580                                                                                                                                                                                                 |
| Gene Symbol                 | SOX13                                                                                                                                                                                                |
| Species                     | Human                                                                                                                                                                                                |
| Immunogen                   | Recombinant protein encompassing a sequence within the center region of human SOX13. The exact sequence is proprietary.                                                                              |
| Product Application Details |                                                                                                                                                                                                      |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation                                                                                                                           |
| Recommended Dilutions       | Western Blot 1:500-1:3000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunoprecipitation 1:100-1:500                                                                                     |




## Images


Western Blot: SOX13 Antibody [NBP1-32847] - (30 ug of whole cell lysate) A: Molt-4 7.5% SDS PAGE diluted at 1:1000



Immunocytochemistry/Immunofluorescence: SOX13 Antibody [NBP1-32847] - Paraformaldehyde-fixed HeLa, using antibody at 1:200 dilution.



Immunoprecipitation: SOX13 Antibody [NBP1-32847] - Immunoprecipitation of SOX13 protein from NT2D1 whole cell extracts using 5 ug of SOX13 antibody [N1C3]. Western blot analysis was performed using SOX13 antibody [N1C3] diluted at 1:500. EasyBlot anti-Rabbit IgG was used as a secondary reagent.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-32847**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00009580-Q01-10ug | Recombinant Human SOX13 GST (N-Term) Protein        |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-32847](http://www.novusbio.com/reviews/submit/NBP1-32847)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

