

Product Datasheet

hnRNP F Antibody - BSA Free

NBP1-33607

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-33607

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-33607



NBP1-33607

hnRNP F Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol
Target Molecular Weight	46 kDa

Product Description


Description	Novus Biologicals Rabbit hnRNP F Antibody - BSA Free (NBP1-33607) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and IP. Anti-hnRNP F Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3185
Gene Symbol	HNRNPF
Species	Human, Mouse
Immunogen	Recombinant protein encompassing a sequence within the center region of human hnRNP F. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Knockdown Validated
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunoprecipitation 1:100-1:500, Immunohistochemistry-Paraffin 1:100-1:1000, Knockdown Validated Validation (PMID: 31814907).

Images


Western Blot: hnRNP F Antibody [NBP1-33607] - Whole cell extract (30 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with hnRNP F antibody [N1N3] diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Immunocytochemistry/Immunofluorescence: hnRNP F Antibody [NBP1-33607] - Paraformaldehyde-fixed A431, using hnRNP F (NBP1-33607) antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: hnRNP F Antibody [NBP1-33607] - Mouse uterus. hnRNP F antibody [N1N3] diluted at 1:500. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Western Blot: hnRNP F Antibody [NBP1-33607] - HepG2 whole cell and nuclear extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with hnRNP F antibody [N1N3] diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Immunohistochemistry-Paraffin: hnRNP F Antibody [NBP1-33607] - A549 Xenograft, using hnRNP F antibody at 1:100 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.



Immunoprecipitation: hnRNP F Antibody [NBP1-33607] - Immunoprecipitation of hnRNP F protein from 293T whole cell extracts using 5 ug of hnRNP F antibody [N1N3]. Western blot analysis was performed using hnRNP F antibody [N1N3]. EasyBlot anti-Rabbit IgG was used as a secondary reagent.



Publications

Xiong L, Gong Y, Liu H et Al. circGlis3 promotes β -cell dysfunction by binding to heterogeneous nuclear ribonucleoprotein F and encoding Glis3-348aa protein iScience 2023-12-09 [PMID: 38226164]

Zhao Q, Liu J, Deng H et al. Targeting Mitochondria-Located circRNA SCAR Alleviates NASH via Reducing mROS Output Cell 2020-09-05 [PMID: 32931733] (WB, Human, Mouse)

Li F, Su M, Zhao H HnRNP-F promotes cell proliferation by regulating TPX2 in bladder cancer Am J Transl Res 2019-11-15 [PMID: 31814907] (KD, IP, WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-33607

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00003185-P01-10ug	Recombinant Human hnRNP F GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-33607

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

