

# Product Datasheet

## MMP-11 Antibody (SL3.05)

### NBP1-40030

Unit Size: 0.05 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-40030](http://www.novusbio.com/NBP1-40030)

Updated 6/20/2021 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-40030](http://www.novusbio.com/reviews/destination/NBP1-40030)



**NBP1-40030**

MMP-11 Antibody (SL3.05)

**Product Information**

|                      |                                                                                        |
|----------------------|----------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.05 ml                                                                                |
| <b>Concentration</b> | 0.2 mg/ml                                                                              |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| <b>Clonality</b>     | Monoclonal                                                                             |
| <b>Clone</b>         | SL3.05                                                                                 |
| <b>Preservative</b>  | 0.09% Sodium Azide                                                                     |
| <b>Isotype</b>       | IgG1 Kappa                                                                             |
| <b>Purity</b>        | Protein G purified                                                                     |
| <b>Buffer</b>        | PBS (pH 7.4) and 0.2% BSA                                                              |

**Product Description**


|                                |                                                                                                                                                                                    |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>                    | Mouse                                                                                                                                                                              |
| <b>Gene ID</b>                 | 4320                                                                                                                                                                               |
| <b>Gene Symbol</b>             | MMP11                                                                                                                                                                              |
| <b>Species</b>                 | Human                                                                                                                                                                              |
| <b>Specificity/Sensitivity</b> | Recognizes human MMP-11 (Stromelysin 3). Preferred for sensitive detection of MMP-11 by Western Blot as well as for staining of routine formalin-fixed, paraffin-embedded tissues. |
| <b>Immunogen</b>               | Recombinant protein corresponding to active form of human MMP-11.                                                                                                                  |

**Product Application Details**

|                              |                                                                                                                      |
|------------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation                               |
| <b>Recommended Dilutions</b> | Western Blot 1:100-1:2000, Immunohistochemistry, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Paraffin 1:200 |

**Images**

Immunohistochemistry-Paraffin: MMP-11 Antibody (SL3.05) [NBP1-40030] - Analysis of anti-MMP11 antibody with human small intestine at dilution 1:200.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-40030](http://www.novusbio.com/reviews/submit/NBP1-40030)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

