

Product Datasheet

Cbx2 Antibody - BSA Free **NBP1-40342**

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-40342

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-40342



NBP1-40342

Cbx2 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)


Product Description	
Host	Rabbit
Gene ID	84733
Gene Symbol	CBX2
Species	Human
Immunogen	The immunogen for this product maps to a region between residue 482 and 532 of human chromobox homolog 2 using the numbering given in entry NP_005180.1 (GeneID 84733).

Product Application Details	
Applications	Immunoprecipitation, Chromatin Immunoprecipitation Sequencing, Western Blot (Negative)
Recommended Dilutions	Immunoprecipitation 2-10 ug/mg lysate, Western Blot (Negative), Chromatin Immunoprecipitation Sequencing 1 ug
Application Notes	Not recommended for Western Blot.




Images

Chromatin Immunoprecipitation Sequencing: Cbx2 Antibody [NBP1-40342] - Localization of CBX2 Binding Sites by ChIP-seq. Chromatin from K562 cells was immunoprecipitated with anti-CBX2 antibody NBP1-40342 and analyzed by DNA sequencing. The figure illustrates the peak distribution of CBX2 binding within a 1.5 Mb region of chromosome 1 as detected using anti-CBX2 antibody NBP1-40342. ChIP-seq validation performed by Diogenode, Denville, NJ.



Immunoprecipitation: Cbx2 Antibody [NBP1-40342] - Detection of human CBX2 by western blot of immunoprecipitates. Samples: Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from HEK293T cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-CBX2 antibody NBP1-40342 (lot NBP1-40342-3) used for IP at 6 ug per reaction. CBX2 was also immunoprecipitated by previous lots of this antibody (lots NBP1-40342-1 and NBP1-40342-2). For blotting immunoprecipitated CBX2, NBP1-40342 was used at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-40342

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00084733-Q01-10ug	Recombinant Human Cbx2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-40342

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

