

# Product Datasheet

## ILT4/CD85d/LILRB2 Antibody (42D1) [PE] NBP1-43772-25tests

Unit Size: 25 Tests

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-43772](http://www.novusbio.com/NBP1-43772)

Updated 7/14/2021 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-43772](http://www.novusbio.com/reviews/destination/NBP1-43772)



**NBP1-43772-25tests**

ILT4/CD85d/LILRB2 Antibody (42D1) [PE]

| Product Information |  |
|---------------------|--|
| Unit Size           | 25 Tests   |
| Concentration       | 0.2 mg/ml  |
| Storage             | Store at 4C. Do not freeze.  |
| Clonality           | Monoclonal   |
| Clone               | 42D1   |
| Preservative        | 0.09% Sodium Azide   |
| Isotype             | IgG2a  |
| Conjugate           | PE   |
| Purity              | Protein A or G purified  |
| Buffer              | PBS [pH 7.2] containing BSA and/or gelatin for protein stabilization |

| Product Description |   |
|---------------------|---|
| Host                | Rat                                       |
| Gene ID             | 10288                                     |
| Gene Symbol         | LILRB2                                    |
| Species             | Human                                     |
| Reactivity Notes    | This antibody is reactive to Human.       |
| Immunogen           | The immunogen for this antibody was ILT4. |

| Product Application Details |   |
|-----------------------------|---|
| Applications                | Flow Cytometry  |
| Recommended Dilutions       | Flow Cytometry  |
| Application Notes           | This 42D1 antibody has been pre-titrated and tested by flow cytometric analysis of normal human peripheral blood cells. This can be used at 20 uL (1 ug)/per test. Cell number should be determined empirically but can range from 10 <sup>5</sup> to 10 <sup>8</sup> cells/test. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-43772](http://www.novusbio.com/reviews/submit/NBP1-43772)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

