

Product Datasheet

DYDC1 Antibody NBP1-47218

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-47218

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-47218



NBP1-47218**DYDC1 Antibody****Product Information**

Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA

Product Description


Host	Goat
Gene ID	143241
Gene Symbol	DYDC1
Species	Human
Immunogen	Peptide with sequence C-RNEDILHSEEATLD corresponding to internal region according to NP_620167.1.

Product Application Details

Applications	Western Blot, Peptide ELISA
Recommended Dilutions	Western Blot 0.5 - 1 ug/ml, Peptide ELISA Detection limit 1:16000
Application Notes	WB: In transfected HEK293 transiently expressing human DYDC1 (myc and DYKDDDDK tagged), a band of approx. 28 kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-myc antibody.

Images

Western Blot: DYDC1 Antibody [NBP1-47218] - Lane A - NBP1-47218 (1ug/ml) staining of HEK293 overexpressing Human DYDC1 lysate (10ug protein in RIPA buffer) Lane B - NBP1-47218 (1ug/ml) staining of HEK293 mock-transfected lysate (10ug protein in RIPA buffer). Lane C - EB00106 anti-MYC Tag (1/1000) staining HEK293 overexpressing Human DYDC1 lysate (10ug protein in RIPA buffer). Primary incubations were for 1 hour. Detected by chemiluminescence.

**Publications**

Deloukas P, Earthrowl ME, Grafham DV et al. The DNA sequence and comparative analysis of human chromosome 10. Nature 2004-05-01 [PMID: 15164054]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-47218

HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control
H00143241-P01-10ug	Recombinant Human DYDC1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-47218

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



