

Product Datasheet

MRG15 Antibody - BSA Free

NBP1-57832

Unit Size: 100ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-57832

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-57832



NBP1-57832

MRG15 Antibody - BSA Free

Product Information

Unit Size	100ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	PBS, 2% Sucrose

Product Description

Description	The addition of 50% glycerol is optional for those storing this antibody at -20C and not aliquoting smaller units. However, please note that glycerol may interrupt some downstream antibody applications and should be added with caution.
Host	Rabbit
Gene ID	10933
Gene Symbol	MORF4L1
Species	Human
Immunogen	Synthetic peptides corresponding to MORF4L1 (mortality factor 4 like 1) The peptide sequence was selected from the middle region of MORF4L1. Peptide sequence YAVNEVVAGIKEYFNVMLGTQLLYKFERPQYAEILADHPDAPMSQ. The peptide sequence for this immunogen was taken from within the described region.

Product Application Details

Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1.0 ug/ml, Immunohistochemistry 1:10-1:500, Immunohistochemistry-Paraffin 1:10-1:500




Images


Western Blot: MRG15 Antibody [NBP1-57832] - HepG2 cell lysate,
Antibody Titration: 5.0-10.0ug/ml




Immunohistochemistry-Paraffin: MRG15 Antibody [NBP1-57832] -
Human urinary bladder tissue, Ganglionic cells (Indicated with Arrows) 4-
8ug/ml.



Immunohistochemistry: MRG15 Antibody [NBP1-57832] - Human cardiac
cell Cellular data: Epithelial cells of renal tubule Antibody Concentration:
4.0-8.0 ug/ml Magnification: 400X.



Immunohistochemistry-Paraffin: MRG15 Antibody [NBP1-57832] -
Human Heart Tissue, antibody concentration 4-8ug/ml. Cells with
positive label: Myocardial cells (indicated with arrows) 400X
magnification.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-57832

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-99031-100ug	Recombinant Human MRG15 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-57832

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

