

Product Datasheet

PMCA4 Antibody NBP1-59481

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-59481

Updated 8/21/2017 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-59481



NBP1-59481

PMCA4 Antibody


| Product Information | |
|---------------------|---|
| Unit Size | 100 ul |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at -20C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Immunogen affinity purified |
| Buffer | PBS and 2% Sucrose |

| Product Description | |
|---------------------|---|
| Host | Rabbit |
| Gene ID | 493 |
| Gene Symbol | ATP2B4 |
| Species | Human |
| Immunogen | Synthetic peptides corresponding to ATP2B4(ATPase, Ca ⁺⁺ transporting, plasma membrane 4) The peptide sequence was selected from the middle region of ATP2B4. Peptide sequence FAGEKFFDIDSGRKAPLHSPPSQHYTIVFNTFVLMQLFNEINSRKIHGEK. |
| Notes | The addition of 50% glycerol is optional for those storing this antibody at -20C and not aliquoting smaller units. However, please note that glycerol may interrupt some downstream antibody applications and should be added with caution. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot |
| Recommended Dilutions | Western Blot 0.2-1 ug/ml |
| Application Notes | This is a rabbit polyclonal antibody against ATP2B4 and was validated on Western blot. |

Images

Western Blot: PMCA4 Antibody [NBP1-59481] - 721_B cell lysate, concentration 0.2-1 ug/ml.





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP1-59481

| | |
|--------------------|--|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP (Horseradish Peroxidase)] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| H00000493-Q01-10ug | Recombinant Human PMCA4 Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-59481

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

