

Product Datasheet

LAIR1 Antibody NBP1-68845

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-68845

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-68845



NBP1-68845**LAIR1 Antibody**

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA


Product Description	
Host	Goat
Gene ID	3903
Gene Symbol	LAIR1
Species	Human
Specificity/Sensitivity	This antibody is expected to recognize both isoforms (NP_002278.1; NP_068352.1).
Immunogen	Peptide with sequence C-NDTEDVSQASPSE corresponding to internal region according to NP_002278.1, NP_068352.1.

Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Peptide ELISA
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry 1-3ug/ml, Immunohistochemistry-Paraffin, Peptide ELISA Detection limit 1:1000
Application Notes	Western blot: Approx 28kDa band observed in Human Lymph Node lysates (calculated MW of 29.8kDa according to NP_068352.1). Primary incubation was 1 hour IHC: Paraffin embedded Human Tonsil.




Images

Western Blot: LAIR1 Antibody [NBP1-68845] - Analysis of LAIR1 in Human Lymph Node lysate (35ug protein in RIPA buffer) using NBP1-68845 at 2ug/ml. Primary incubation was 1 hour. Detected by chemiluminescence.



Immunohistochemistry-Paraffin: LAIR1 Antibody [NBP1-68845] - 3.75ug/ml staining of paraffin embedded Human Tonsil. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-68845

HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control
NBP1-84590PEP	LAIR1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-68845

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



