

Product Datasheet

IKK epsilon/IKBKE Antibody - BSA Free NBP1-76800

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-76800

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-76800



NBP1-76800

IKK epsilon/IKBKE Antibody - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	80 kDa

Product Description


Host	Rabbit
Gene ID	9641
Gene Symbol	IKBKE
Species	Human, Mouse
Specificity/Sensitivity	It has no cross response to IKKa, IKKb, or IKKy. Human IKK epsilon/IKBKE has 2 isoforms, including isoform 1 (716aa, 80kD) and isoform 2 (631aa, 70kD). Mouse IKK epsilon/IKBKE has one isoform (717aa, 81kD) and Rat IKK epsilon/IKBKE also has one isoform (717aa, 80kD). NBP1-76800 can detect human, mouse and possibly rat.
Immunogen	Antibody was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of human IKK epsilon/IKBKE. The immunogen is located within the last 50 amino acids of IKK epsilon/IKBKE. Amino Acid Sequence: NRIIERLNRVPAPPDV

Product Application Details


Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1-2 ug/ml, ELISA 1:100-1:2000, Immunohistochemistry 10 ug/ml, Immunocytochemistry/ Immunofluorescence 20 ug/ml, Immunohistochemistry-Paraffin 10 ug/ml

Images


Western Blot: IKK epsilon/IKBKE Antibody [NBP1-76800] - Western Blot Validation in Jurkat Cell Lysate. Loading: 15 ug of lysates per lane. Antibodies: IKK epsilon/IKBKE NBP1-76800 (2 ug/mL), 1h incubation at RT in 5% NFDN/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.




Immunocytochemistry/Immunofluorescence: IKK epsilon/IKBKE Antibody [NBP1-76800] - Hela cells with IKK epsilon antibody at 20 ug/mL. Green: IKK epsilon Antibody (NBP1-76800). Blue: DAPI staining.




Immunohistochemistry-Paraffin: IKK epsilon/IKBKE Antibody [NBP1-76800] - Human pancreas tissue with IKBKE antibody at 10 ug/ml.




Western Blot: IKK epsilon/IKBKE Antibody [NBP1-76800] - Induced Expression Validation in Mouse Cell Line. Loading: 15 ug of RAW 264.7 cell lysates per lane. Antibodies: IKK epsilon/IKBKE NBP1-76800 (1 ug/ml), 1h incubation at RT in 5% NFDN/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution. Cells were treated with LPS (0.3 ug/ml) for 3 hrs (lane 2) and 6 hrs (lane 3) or with control (lane 1).



Western Blot: IKK epsilon/IKBKE Antibody [NBP1-76800] - Western Blot Validation in Human and Mouse Cell Lines. Loading: 15 ug of lysates per lane. Antibodies: IKK epsilon/IKBKE NBP1-76800 (1 ug/mL), 1h incubation at RT in 5% NFDN/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.



Western Blot: IKK epsilon/IKBKE Antibody [NBP1-76800] - Western Blot Validation in Mouse Tissue Lysates. Loading: 15 ug of lysates per lane. Antibodies: IKK epsilon/IKBKE NBP1-76800 (2 ug/mL), 1h incubation at RT in 5% NFDm/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.



Publications

Peant B, Diallo JS, Lessard L et al. Regulation of IkappaB kinase epsilon expression by the androgen receptor and the nuclear factor-kappaB transcription factor in prostate cancer. *Mol Cancer Res*. 2007-01-01 [PMID: 17259348] (WB, Human)

Details:

Imgenex antibodies cited [WB (human androgen-dependent prostate cancer cell lines), Figs. 1A, 3B, 4, 5A,B]: 1. IMG-5571 (IKK-epsilon) 2. IMG-139A (TBK-1).



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-76800

NBP1-76800PEP	IKK epsilon/IKBKE Antibody Blocking Peptide
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-76800

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

