

Product Datasheet

DFF40/CAD Antibody - BSA Free NBP1-76884

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-76884

Updated 4/24/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-76884



NBP1-76884


DFF40/CAD Antibody - BSA Free

| Product Information | |
|-----------------------------|--|
| Unit Size | 0.1 mg |
| Concentration | 1 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.02% Sodium Azide |
| Isotype | IgG |
| Purity | Immunogen affinity purified |
| Buffer | PBS |
| Target Molecular Weight | 40 kDa |
| Product Description | |
| Host | Rabbit |
| Gene ID | 1677 |
| Gene Symbol | DFFB |
| Species | Human |
| Immunogen | Antibody was raised against a 16 amino acid peptide near the amino terminus of human DFF40. The immunogen is located within the first 50 amino acids of DFF40. |
| Product Application Details | |
| Applications | Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence |
| Recommended Dilutions | Western Blot 0.5 ug/mL, ELISA 1:100-1:2000, Immunocytochemistry/ Immunofluorescence 20 ug/mL |
| Application Notes | DFF40 antibody can be used for detection of DFF40 by Western blot at 1:500 to 1:2000 dilution. A 40 kDa band can be detected. |




Images


Western Blot: DFF40/CAD Antibody [NBP1-76884] - Jurkat (J) and K562 (K) whole cell lysate with DFF40 antibody at 1:500 dilution.



Immunocytochemistry/Immunofluorescence: DFF40/CAD Antibody [NBP1-76884] - Jurkat cells with DFF40 antibody at 5 ug/mL.



Immunocytochemistry/Immunofluorescence: DFF40/CAD Antibody [NBP1-76884] - Immunofluorescence of DFF40 in Jurkat cells with DFF40 antibody at 20 ug/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-76884

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

