

Product Datasheet

Phosphoserine Antibody [Biotin] NBP1-77727

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-77727

Updated 9/9/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-77727



NBP1-77727**Phosphoserine Antibody [Biotin]**

Product Information	
Unit Size	0.1 mg
Concentration	0.25 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Conjugate	Biotin
Purity	Immunogen affinity purified
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 40% (v/v) Glycerol
Product Description	
Host	Rabbit
Species	Non-species specific
Specificity/Sensitivity	The antibody reacts with free phosphoserine, phosphoserine conjugated to carriers such as thyroglobulin or BSA, and detects the presence of phosphoserine in proteins of both unstimulated and stimulated cell lysates. It shows minimal reactivity by ELISA and western blotting with proteins containing phosphotyrosine or phosphothreonine
Immunogen	Phosphoserine Antibody was prepared from whole rabbit serum produced by repeated immunizations with phosphoserine conjugated KLH.
Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500, ELISA 1:5000, Immunohistochemistry 1:10-1:500, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Paraffin 1:10-1:500
Application Notes	This affinity-purified antibody is suitable for use in ELISA and western blotting assays. The antibody reacts specifically with phosphoserine and shows minimal reactivity by ELISA and western blotting with proteins containing phosphotyrosine or phosphothreonine residues. The antibody reacts with free phosphoserine, phospho-serine conjugated - carriers such as thyroglobulin or BSA, and detects the presence of phosphoserine in proteins of both unstimulated and stimulated cell lysates. Although not tested, this antibody is likely functional in RIA, flow cytometry and immunohistochemistry.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-77727

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

