

Product Datasheet

RFX5 Antibody NBP1-77821

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-77821

Updated 5/3/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-77821



NBP1-77821**RFX5 Antibody**

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Sodium Azide
Isotype	IgG
Purity	Ion exchange chromatography
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Product Description	
Description	This antibody was prepared from monospecific antiserum by delipidation, salt fractionation and ion exchange chromatography.
Host	Rabbit
Gene ID	5993
Gene Symbol	RFX5
Species	Human
Specificity/Sensitivity	Anti-RFX5 (N-terminal specific) may react with unknown minor bands at 55 kDa and 85 kDa.
Immunogen	RFX5 (N-terminal specific) peptide corresponding to a region near the N-terminus of the human protein conjugated to Keyhole Limpet Hemocyanin (KLH).

Product Application Details	
Applications	Western Blot, ELISA, Gel Super Shift Assays
Recommended Dilutions	Western Blot 1:500-1:3000, ELISA 1:5000-1:25000, Gel Super Shift Assays
Application Notes	This product was assayed by immunoblot and found to be reactive against RFX5 from a variety of fibroblast and B-cell lysates at a dilution of 1:500 followed by reaction with Peroxidase conjugated Affinity Purified anti-Rabbit IgG [H&L] Anti-RFX5 is suitable for the detection by immunoblot of human RFX5 a shows a 75 kDa band. This product was also tested in a gel supershift assay and found to be reactive against RFX5 complexes using 0.5 - 1.0 ul per assay. This product was assayed by immunoblot and found to be reactive against RFX5 from a variety of fibroblast and B-cell lysates at a dilution of 1:500 followed by reaction with Peroxidase conjugated Affinity Purified anti-Rabbit IgG [H&L] Anti-RFX5 is suitable for the detection by immunoblot of human RFX5 a shows a 75 kDa band. This product was also tested in a gel supershift assay and found to be reactive against RFX5 complexes using 0.5 - 1.0 ul per assay.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-77821

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

