

# Product Datasheet

## GSAP Antibody [FITC] NBP1-78400F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-78400F](http://www.novusbio.com/NBP1-78400F)

Updated 10/29/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-78400F](http://www.novusbio.com/reviews/destination/NBP1-78400F)



**NBP1-78400F**

GSAP Antibody [FITC]

| Product Information         |                                                                                                                       |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                           |
| Storage                     | Store at 4C in the dark.                                                                                              |
| Clonality                   | Polyclonal                                                                                                            |
| Preservative                | 0.05% Sodium Azide                                                                                                    |
| Isotype                     | IgG                                                                                                                   |
| Conjugate                   | FITC                                                                                                                  |
| Purity                      | Immunogen affinity purified                                                                                           |
| Buffer                      | PBS                                                                                                                   |
| Product Description         |                                                                                                                       |
| Host                        | Rabbit                                                                                                                |
| Gene ID                     | 54103                                                                                                                 |
| Gene Symbol                 | PION                                                                                                                  |
| Species                     | Human                                                                                                                 |
| Reactivity Notes            | Human.                                                                                                                |
| Immunogen                   | A synthetic peptide made to an internal portion of the human GSAP protein (between residues 150-200) [UniProt A4D1B5] |
| Product Application Details |                                                                                                                       |
| Applications                | Western Blot                                                                                                          |
| Recommended Dilutions       | Western Blot                                                                                                          |
| Application Notes           | This GSAP antibody is useful for Western Blot where a band is observed ~65 kDa representing Isoform 2/3 of GSAP.      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-78400F**

---

|               |                                   |
|---------------|-----------------------------------|
| NBP2-24892    | Rabbit IgG Isotype Control [FITC] |
| NBP1-78400PEP | GSAP Antibody Blocking Peptide    |
| 664-LI-025    | LIGHT/TNFSF14 [Unconjugated]      |
| NBP1-78376PEP | GSAP Antibody Blocking Peptide    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-78400F](http://www.novusbio.com/reviews/submit/NBP1-78400F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

