

# Product Datasheet

## HNRPA3 Antibody - BSA Free

### NBP1-80486

Unit Size: 100ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-80486](http://www.novusbio.com/NBP1-80486)

Updated 2/21/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-80486](http://www.novusbio.com/reviews/destination/NBP1-80486)



**NBP1-80486**

HNRPA3 Antibody - BSA Free

**Product Information**

|                      |                                                                                        |
|----------------------|----------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100ul                                                                                  |
| <b>Concentration</b> | 1 mg/ml                                                                                |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| <b>Clonality</b>     | Polyclonal                                                                             |
| <b>Preservative</b>  | 0.09% Sodium Azide                                                                     |
| <b>Isotype</b>       | IgG                                                                                    |
| <b>Purity</b>        | Protein A purified                                                                     |
| <b>Buffer</b>        | PBS, 2% Sucrose                                                                        |

**Product Description**


|                    |                                                                                                                                                                                                                                             |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | The addition of 50% glycerol is optional for those storing this antibody at -20C and not aliquoting smaller units. However, please note that glycerol may interrupt some downstream antibody applications and should be added with caution. |
| <b>Host</b>        | Rabbit                                                                                                                                                                                                                                      |
| <b>Gene ID</b>     | 220988                                                                                                                                                                                                                                      |
| <b>Gene Symbol</b> | HNRNPA3                                                                                                                                                                                                                                     |
| <b>Species</b>     | Human                                                                                                                                                                                                                                       |
| <b>Immunogen</b>   | Synthetic peptide directed towards the N terminal of human HNRPA3. Peptide sequence<br>MEVKPPPGRPQPDSGRRRRRRGEEGHDPKEPEQLRKLFIGGLSFETTTDS.<br>The peptide sequence for this immunogen was taken from within the described region.           |

**Product Application Details**


|                              |                                                                                                   |
|------------------------------|---------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin                                 |
| <b>Recommended Dilutions</b> | Western Blot 1.0 ug/ml, Immunohistochemistry 1:10-1:500, Immunohistochemistry-Paraffin 1:10-1:500 |

**Images**


Western Blot: HNRPA3 Antibody [NBP1-80486] - Host: Rabbit. Target: HNRNPA3. Positive control (+): HeLa Cell Lysate (HL). Negative control (-): Human Fetal Heart (HE). Antibody concentration: 3ug/ml




Immunohistochemistry: HNRPA3 Antibody [NBP1-80486] - Human Liver  
 Cellular data: Hepatocytes Antibody Concentration: 4.0-8.0 ug/ml  
 Magnification: 400X.




Western Blot: HNRPA3 Antibody [NBP1-80486] - Jurkat cell lysate,  
 Antibody Titration: 1.25ug/ml




Western Blot: HNRPA3 Antibody [NBP1-80486] - Sample Tissue: Jurkat,  
 Lane A: Primary Antibody, Lane B: Primary Antibody + Blocking Peptide,  
 Primary Antibody Concentration: 0.625ug/mL, Peptide Concentration:  
 1.0ug/mL, Lysate Quantity: 25ug/lane, Gel Concentration: 12%



Western Blot: HNRPA3 Antibody [NBP1-80486] - Antibody Titration: 1  
 ug/ml Human 721\_B.




Western Blot: HNRPA3 Antibody [NBP1-80486] - Antibody Titration: 1 ug/ml Human Daudi.




Western Blot: HNRPA3 Antibody [NBP1-80486] - Sample Tissue: Human Jurkat Antibody Dilution: 1.0 ug/ml

### HNRPA3



Immunohistochemistry-Paraffin: HNRPA3 Antibody [NBP1-80486] - Human Liver Tissue, antibody concentration 4-8ug/ml. Cells with positive label: Hepatocytes (indicated with arrows) 400X magnification.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-80486**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP2-56989PEP | HNRPA3 Recombinant Protein Antigen                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-80486](http://www.novusbio.com/reviews/submit/NBP1-80486)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

