

Product Datasheet

TSG101 Antibody - BSA Free

NBP1-80659

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 2 Publications: 11

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-80659

Updated 4/13/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-80659



NBP1-80659

TSG101 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Rabbit
Gene ID	7251
Gene Symbol	TSG101
Species	Human, Mouse, Rat
Marker	Exosome Marker
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: RMKEEMDRAQAELNALKRTEEDLKKGHQKLEEMVTRLDQEVAEVDKNIELLKK KDEELSSALEKMENQSENNDIDEVIIPTAPLYKQILNLYAEENAIEDTIFYLGEALR RGVIDLDDVFLKHVRLLS


Product Application Details	
Applications	Western Blot, Simple Western, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockdown Validated
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Simple Western 1:50, Immunohistochemistry 1:200 - 1:500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:200 - 1:500, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization, Use PFA/Triton X-100.

Images


Western Blot: TSG101 Antibody [NBP1-80659] - Analysis in MCF-7 cells transfected with control siRNA, target specific siRNA probe #1 and #2, using Anti-TSG101 antibody. Remaining relative intensity is presented. Loading control: Anti-GAPDH.




Immunocytochemistry/Immunofluorescence: TSG101 Antibody [NBP1-80659] - Staining of human cell line A-431 shows localization to plasma membrane & cytosol. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: TSG101 Antibody [NBP1-80659] - Staining of human fallopian tube shows strong cytoplasmic positivity in glandular cells.




Western Blot: TSG101 Antibody [NBP1-80659] - Analysis in mouse cell line NIH-3T3 and rat cell line NBT-II.




Immunohistochemistry-Paraffin: TSG101 Antibody [NBP1-80659] - Staining of human prostate shows moderate cytoplasmic positivity in glandular and in smooth muscle cells.



Immunohistochemistry-Paraffin: TSG101 Antibody [NBP1-80659] - Staining of human cerebral cortex shows moderate cytoplasmic positivity in neurons.



Immunohistochemistry-Paraffin: TSG101 Antibody [NBP1-80659] - Staining of human small intestine shows strong cytoplasmic positivity in glandular cells.



Publications

Jiao W, Hao J, Liu JM et al. Mesenchymal stem cells-derived extracellular vesicle-incorporated H19 attenuates cardiac remodeling in rats with heart failure The Kaohsiung journal of medical sciences 2023-10-27 [PMID: 37885317] (Western Blot, Rat)

Zar M, Campodonico J, Cosentino N et al. Plasma Exosome Profile in ST-Elevation Myocardial Infarction Patients with and without Out-of-Hospital Cardiac Arrest International Journal of Molecular Sciences 2021-07-28 [PMID: 34360827] (Western Blot)

Majewska A, Brodaczewska K, Filipiak-Duliban A et al. miRNA Pattern in Hypoxic Microenvironment of Kidney Cancer-Role of PTEN Biomolecules 2022-05-11 [PMID: 35625614] (Western Blot)

Carotenuto P, Romano A, Barbato A et al. Targeting the MITF/APAF-1 axis as salvage therapy for MAPK inhibitors in resistant melanoma Cell reports 2022-11-08 [PMID: 36351409] (WB, ICC/IF, Human)

Kim EJ, Kim JY, Kim SO et al. The oncogenic JAG1 intracellular domain is a transcriptional cofactor that acts in concert with DDX17/SMAD3/TGIF2 Cell reports 2022-11-22 [PMID: 36417870] (WB, Human)

Hua Y, Dai C, He Q Et al. Autoantibody panel on small extracellular vesicles for the early detection of lung cancer Clin Immunol 2022-11-10 [PMID: 36356847] (WB, Human)

Details:

Citation using the PE version of this antibody.

Pauwels MJ, Xie J, Ceroi A et al. Choroid plexus-derived extracellular vesicles exhibit brain targeting characteristics Biomaterials 2022-11-01 [PMID: 36302306] (WB, Mouse)

Vandendriessche C, Balusu S, Van Cauwenberghe C Et al. Importance of extracellular vesicle secretion at the blood-cerebrospinal fluid interface in the pathogenesis of Alzheimer's disease Acta neuropathologica communications 2021-08-23 [PMID: 34425919] (WB, Mouse)

Millan C, Prause L, Vallmajo-Martin Q et al. Extracellular Vesicles from 3D Engineered Microtissues Harbor Disease-Related Cargo Absent in EVs from 2D Cultures Advanced healthcare materials 2021-04-22 [PMID: 33890421]

Pospichalova V, Svoboda J, Dave Z et al. Simplified protocol for flow cytometry analysis of fluorescently labeled exosomes and microvesicles using dedicated flow cytometer. J Extracell Vesicles 2015-01-01 [PMID: 25833224] (WB)

Gross JC, Chaudhary V, Bartscherer K, Boutros M. Active Wnt proteins are secreted on exosomes. Nat Cell Biol 2012-10-01 [PMID: 22983114]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-80659

NBP1-80659PEP	TSG101 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-80659

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

