

# Product Datasheet

## **B3GAT2 Antibody** **NBP1-81828**

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-81828](http://www.novusbio.com/NBP1-81828)

Updated 2/12/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-81828](http://www.novusbio.com/reviews/destination/NBP1-81828)



**NBP1-81828****B3GAT2 Antibody****Product Information**

|                      |                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                                           |
| <b>Concentration</b> | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| <b>Clonality</b>     | Polyclonal                                                                                                       |
| <b>Preservative</b>  | 0.02% Sodium Azide                                                                                               |
| <b>Isotype</b>       | IgG                                                                                                              |
| <b>Purity</b>        | Immunogen affinity purified                                                                                      |
| <b>Buffer</b>        | PBS (pH 7.2) and 40% Glycerol                                                                                    |

**Product Description**


|                    |                                                                                                                                                                                                                    |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>        | Rabbit                                                                                                                                                                                                             |
| <b>Gene ID</b>     | 135152                                                                                                                                                                                                             |
| <b>Gene Symbol</b> | B3GAT2                                                                                                                                                                                                             |
| <b>Species</b>     | Human                                                                                                                                                                                                              |
| <b>Immunogen</b>   | This antibody was developed against Recombinant Protein corresponding to amino acids:<br>PVGLVGRRYERPLVENGKVVGWYTGWRA DRPFAIDMAGFAVSLQVILSNPK<br>AVFKRRGSQPGMQESDFLKQITTVEELEPKANNCTKVLVWHTRTEKVNLANEP<br>KYHLDTVK |

**Product Application Details**

|                              |                                                                             |
|------------------------------|-----------------------------------------------------------------------------|
| <b>Applications</b>          | Immunohistochemistry, Immunohistochemistry-Paraffin                         |
| <b>Recommended Dilutions</b> | Immunohistochemistry 1:20 - 1:50, Immunohistochemistry-Paraffin 1:20 - 1:50 |
| <b>Application Notes</b>     | For IHC-Paraffin, HIER pH 6 retrieval is recommended.                       |

**Images**

Immunohistochemistry-Paraffin: B3GAT2 Antibody [NBP1-81828] - Staining of human corpus, uterine shows strong cytoplasmic and membranous positivity in glandular cells.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-81828](http://www.novusbio.com/reviews/submit/NBP1-81828)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

