

Product Datasheet

Human Resistin ELISA Kit (Colorimetric) NBP1-83744

Unit Size: 96 Tests

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-83744

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-83744



NBP1-83744**Human Resistin ELISA Kit (Colorimetric)****Product Information**

Unit Size	96 Tests
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.

Product Description


Description	This kit is an enzyme-linked immunosorbent assay for quantitative detection.
Gene ID	56729
Gene Symbol	RETN
Species	Human
Specificity/Sensitivity	The interference of circulating factors of the immune system was evaluated by spiking these proteins at physiologically relevant concentrations into a human Resistin positive serum. There was no crossreactivity detected.
Kit Components	Aluminium pouch(es) with a Microwell Plate coated with polyclonal antibody to human Resistin, Biotin-Conjugate anti-human Resistin polyclonal antibody, Streptavidin-HRP, Human Resistin Standard lyophilized, 4000 pg/ml upon reconstitution, Assay Buffer Concentrate 20x (PBS with 1% Tween 20 and 10% BSA), Wash Buffer Concentrate 20x (PBS with 1% Tween 20), Substrate Solution (tetramethyl-benzidine), Stop Solution (1M Phosphoric acid), Adhesive Films
Standard Curve Range	31 - 2,000 pg/mL
Sensitivity	3.1 pg/mL
Inter Assay	CV% < 8.1%
Intra Assay	CV% < 5.1%
Assay Type	Sandwich ELISA
Suitable Sample Type	Cell culture supernatant, serum, plasma (EDTA, citrate, heparin), buffered solution
Sample Volume	50 uL

Product Application Details

Applications	ELISA
Recommended Dilutions	ELISA

Images

ELISA: Resistin ELISA Kit [NBP1-83744] - Standard curve





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-83744

NBP2-61343-5ug	Recombinant Rat Resistin Protein
210-TA-005	TNF-alpha [Unconjugated]
DRSN00	Resistin [HRP]
M6000B-1	IL-6 [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-83744

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

