

Product Datasheet

SAMM50 Antibody - BSA Free

NBP1-84509

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-84509

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-84509



NBP1-84509

SAMM50 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Rabbit SAMM50 Antibody - BSA Free (NBP1-84509) is a polyclonal antibody validated for use in IHC, WB, ICC/IF, Simple Western and IP. Anti-SAMM50 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	25813
Gene Symbol	SAMM50
Species	Human, Mouse, Rat
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: HARSLEPLPSSGPDFGGLGEEAEFVEVEPEAKQEILENKDVVVQHVHFDGLGR TKDDIICEIGDVFKAKNLIEVMRKSHEAREKLLRLGIFRQVDVLID

Product Application Details	
Applications	Western Blot, Simple Western, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Simple Western 1:4, Immunohistochemistry 1:1000 - 1:2500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunoprecipitation Reported in scientific literature (PMID:26308446)., Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	ICC/IF Fixation Permeabilization: Use PFA/Triton X-100. IHC-Paraffin HIER pH6 retrieval is recommended. See Simple Western Antibody Database for Simple Western validation: Tested in RT-4 and U-251MG, separated by Size, antibody dilution of 1:4, apparent MW was 56 kDa




Images


Immunocytochemistry/Immunofluorescence: SAMM50 Antibody [NBP1-84509] - Staining of human cell line A-431 shows localization to mitochondria. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: SAMM50 Antibody [NBP1-84509] - Staining of human parathyroid gland shows very strong granular cytoplasmic positivity in glandular cells.




Immunohistochemistry-Paraffin: SAMM50 Antibody [NBP1-84509] - Staining of human small intestine shows strong cytoplasmic granular positivity in glandular cells.




Immunohistochemistry-Paraffin: SAMM50 Antibody [NBP1-84509] - Staining of human kidney shows very strong cytoplasmic granular positivity in cells in tubules.




Immunohistochemistry-Paraffin: SAMM50 Antibody [NBP1-84509] - Staining of human pancreas shows moderate cytoplasmic granular positivity in exocrine glandular cells.




Simple Western: SAMM50 Antibody [NBP1-84509] - Simple Western lane view shows a specific band for SAMM50 in 0.2 mg/ml of RT-4 (Left) and U-251MG (Right) lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.




Simple Western: SAMM50 Antibody [NBP1-84509] - Electropherogram image(s) of corresponding Simple Western lane view. SAMM50 antibody was used at 1:4 dilution on RT-4 and U-251MG lysate(s).




Staining of human testis shows strong cytoplasmic granular positivity in Leydig cells.




Analysis in mouse cell line NIH-3T3 and rat cell line NBT-II.



Analysis in human cell line HL-60.




Analysis in human kidney and pancreas tissues using HPA034537 antibody. Corresponding SAMM50 RNA-seq data are presented for the same tissues.



SAMM50 in Kidney

SAMM50 in Pancreas



Publications

Takamitsu E, Otsuka M, Haebara T et al. Identification of Human N-Myristoylated Proteins from Human Complementary DNA Resources by Cell-Free and Cellular Metabolic Labeling Analyses. PLoS One 2015-01-01 [PMID: 26308446] (IP, ICC/IF, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-84509

NBP1-84509PEP	SAMM50 Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-84509

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

