

Product Datasheet

DHX58 Antibody - BSA Free

NBP1-85348

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-85348

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-85348



NBP1-85348

DHX58 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Rabbit
Gene ID	79132
Gene Symbol	DHX58
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: AINHVLQLCANLDTWCIMSPQNCCPQLQEHSQQPCKQYNLCHRRSQDPFGDL LKKLMDQIHDHLEMPELSRKFGTQMYEQQVVKLS

Product Application Details


Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:500 - 1:1000, Immunohistochemistry-Paraffin 1:500 - 1:1000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Western Blot: DHX58 Antibody [NBP1-85348] - (A) HeLa and (B) C2C12 cells were treated with the indicated siRNAs for 2 days and cellular lysates were subjected to Western blot analysis using anti-rapsyn antibody. GAPDH was used as loading control. Rapsyn protein levels are reduced in NUP88-depleted cells. Quantification of the respective expression levels of rapsyn after transfection of HeLa and C2C12 cells with the indicated siRNAs and shRNA-mediated depletion of NUP88 in C2C12 cells. Blots from three independent experiments for each condition were analyzed. Data present mean \pm SEM. P-values: *** <0.001 ; ** <0.01 ; * <0.05 ; t-test, one-tailed. Image collected and cropped by CiteAb from the following publication (<https://dx.plos.org/10.1371/journal.pgen.1007845>) licensed under a CC-BY license.



Immunohistochemistry: DHX58 Antibody [NBP1-85348] - Bright-field images of histological muscle sections from individual B.II.2 and a control fetus stained with anti-rapsyn antibodies (brown). Nuclei were visualized by hematoxylin. Image collected and cropped by CiteAb from the following publication (<https://dx.plos.org/10.1371/journal.pgen.1007845>) licensed under a CC-BY license.




Western Blot: DHX58 Antibody [NBP1-85348] - Lane 1: Marker [kDa] 230, 130, 95, 72, 56, 36, 28, 17, 11. Lane 2: Human cell line RT-163



Western Blot: DHX58 Antibody [NBP1-85348] - Rig-I Like Receptor (RLR) gene expression in stable cell lines encompassing ST-RLRs. Western blot analysis was performed using (D) anti-LGP2 in ST-RLR cells. Cells were transfected with pIFNB-FLuc, pTK-Rluc and either mock, poly(I:C) or 5'3P. Image collected and cropped by CiteAb from the following publication (<https://elifesciences.org/articles/11275>) licensed under a CC-BY license.



Immunohistochemistry-Paraffin: DHX58 Antibody [NBP1-85348] - Staining of human kidney shows strong cytoplasmic positivity in cells in tubules.



Publications

David RYS, Combredet C, Sismeiro O. Comparative analysis of viral RNA signatures on different RIG-I-like receptors. Elife. 2016-03-24 [PMID: 27011352] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-85348

NBP1-85348PEP	DHX58 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-85348

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

