

# Product Datasheet

## Eps8 Antibody - BSA Free NBP1-86150

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-86150](http://www.novusbio.com/NBP1-86150)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-86150](http://www.novusbio.com/reviews/destination/NBP1-86150)



**NBP1-86150**


Eps8 Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                                                          |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                   |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                                                                                         |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                                                                                                                                   |
| Clonality                   | Polyclonal                                                                                                                                                                                                                               |
| Preservative                | 0.02% Sodium Azide                                                                                                                                                                                                                       |
| Isotype                     | IgG                                                                                                                                                                                                                                      |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                                              |
| Buffer                      | PBS (pH 7.2) and 40% Glycerol                                                                                                                                                                                                            |
| Product Description         |                                                                                                                                                                                                                                          |
| Description                 | Novus Biologicals Rabbit Eps8 Antibody - BSA Free (NBP1-86150) is a polyclonal antibody validated for use in IHC and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.                                         |
| Host                        | Rabbit                                                                                                                                                                                                                                   |
| Gene ID                     | 2059                                                                                                                                                                                                                                     |
| Gene Symbol                 | EPS8                                                                                                                                                                                                                                     |
| Species                     | Human                                                                                                                                                                                                                                    |
| Reactivity Notes            | Immunogen displays the following percentage of sequence identity for non-tested species: Mouse (80%), Rat (80%)                                                                                                                          |
| Immunogen                   | This antibody was developed against Recombinant Protein corresponding to amino acids:<br>PKEQFIPPYVPRFRNGWEPPMLNFMGATMEQDLYQLAESVANVAEHQRKQEIKRLSTEHSSVSEYHPADGYAFSSNIYTRGSHLDQGEAAVAFKPTSNRHIDRNYEPLKTQPKKYAKSKYDFVARNNSELSVLKDDILEILDD |
| Product Application Details |                                                                                                                                                                                                                                          |
| Applications                | Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry                                                                                                                                             |
| Recommended Dilutions       | Immunohistochemistry 1:20 - 1:50, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:20 - 1:50                                                                                                        |
| Application Notes           | For IHC-Paraffin, HIER pH 6 retrieval is recommended.<br>Immunocytochemistry/Immunofluorescence Fixation Permeabilization: Use PFA/Triton X-100.                                                                                         |




## Images


Immunocytochemistry/Immunofluorescence: Eps8 Antibody [NBP1-86150] - Immunofluorescent staining of human cell line U-251 MG shows localization to the Golgi apparatus.




Staining of human kidney shows moderate cytoplasmic positivity in cells in tubules.




Staining of human rectum shows moderate cytoplasmic positivity in glandular cells.



Staining of human pancreas shows moderate cytoplasmic positivity in exocrine glandular cells.



Staining of human placenta shows moderate cytoplasmic positivity in trophoblastic cells.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP1-86150**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| NBP1-86150PEP | Eps8 Recombinant Protein Antigen                    |
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-86150](http://www.novusbio.com/reviews/submit/NBP1-86150)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

