

Product Datasheet

CDK5RAP3 Antibody - BSA Free

NBP1-86778

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-86778

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-86778



NBP1-86778

CDK5RAP3 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Rabbit
Gene ID	80279
Gene Symbol	CDK5RAP3
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: GDFGVEAVSEGTDSGISAEAAAGIDWGIFPESDSKDPGGDGIDWGDDAVALQIT VLEAGTQAPEGVARGPDALTLLTYTETRNQFLDELMELEIFLAQRAVELSEED VLSVSQFQ

Product Application Details


Applications	Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry 1:1000 - 1:2500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. Immunocytochemistry/Immunofluorescence Fixation Permeabilization: Use PFA/Triton X-100.

Images


Immunocytochemistry/Immunofluorescence: CDK5RAP3 Antibody [NBP1-86778] - Staining of human cell line U-2 OS shows localization to cytosol.




Immunohistochemistry-Paraffin: CDK5RAP3 Antibody [NBP1-86778] - Staining of human skeletal muscle shows weak cytoplasmic positivity in myocytes.




Immunohistochemistry: CDK5RAP3 Antibody [NBP1-86778] - Staining of human kidney shows strong cytoplasmic positivity in tubule cells.




Immunohistochemistry-Paraffin: CDK5RAP3 Antibody [NBP1-86778] - Staining of human epididymis shows strong cytoplasmic positivity in glandular cells.




Immunohistochemistry-Paraffin: CDK5RAP3 Antibody [NBP1-86778] - Analysis in human epididymis and skeletal muscle tissues using NBP1-86778 antibody. Corresponding CDK5RAP3 RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: CDK5RAP3 Antibody [NBP1-86778] - Staining of human epididymis shows strong cytoplasmic positivity in glandular cells.




Immunohistochemistry-Paraffin: CDK5RAP3 Antibody [NBP1-86778] - Staining of human fallopian tube shows strong cytoplasmic positivity in glandular cells.



Immunohistochemistry-Paraffin: CDK5RAP3 Antibody [NBP1-86778] - Staining of human prostate shows strong cytoplasmic positivity in glandular cells.



Immunohistochemistry-Paraffin: CDK5RAP3 Antibody [NBP1-86778] - Staining of human testis shows strong cytoplasmic positivity in cells in seminiferous ducts and leydig cells.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-86778

NBP1-86778PEP	CDK5RAP3 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-86778

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

