

Product Datasheet

Cdk6 Antibody - BSA Free NBP1-87262

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-87262

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-87262



NBP1-87262

Cdk6 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Description	Novus Biologicals Rabbit Cdk6 Antibody - BSA Free (NBP1-87262) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-Cdk6 Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1021
Gene Symbol	CDK6
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: FRGSSDVVDQLGKILDVIGLPGEEDWPRDVALPRQAFHKSQAQPIEFVTDIDEL GKDLLLKCLTFNPAKRISAYSALSHPYFQDLERCKENLDSHLPPSQNTSELN

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images


Immunohistochemistry-Paraffin: Cdk6 Antibody [NBP1-87262] - Analysis in human placenta and skeletal muscle tissues. Corresponding CDK6 RNA-seq data are presented for the same tissues.




Western Blot: Cdk6 Antibody [NBP1-87262] - Analysis in Rh30 cells transfected with control siRNA, target specific siRNA probe #1 and #2. Remaining relative intensity is presented. Loading control: Anti-PPIB.




Immunocytochemistry/Immunofluorescence: Cdk6 Antibody [NBP1-87262] - Staining of human cell line U-251 MG shows localization to nucleoplasm and cytosol. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: Cdk6 Antibody [NBP1-87262] - Staining of human skeletal muscle shows no positivity in myocytes as expected.




Western Blot: Cdk6 Antibody [NBP1-87262] - Analysis in human cell lines U2OS and A-549. Corresponding RNA-seq data are presented for the same cell lines. Loading control: Anti-COX4I1.




Immunohistochemistry-Paraffin: Cdk6 Antibody [NBP1-87262] - Staining of human lymph node shows moderate to strong cytoplasmic positivity in germinal center cells.



Immunohistochemistry-Paraffin: Cdk6 Antibody [NBP1-87262] - Staining of human testis shows moderate nuclear positivity in cells in seminiferous ducts.



Immunohistochemistry-Paraffin: Cdk6 Antibody [NBP1-87262] - Staining of human placenta shows strong nuclear positivity in trophoblastic cells.



Publications

Sikorski K, Mehta A et al. A high-throughput pipeline for validation of antibodies. Nat Methods 2018-01-11 [PMID: 30377371] (Human)

Details:

Antibody validation based on denaturing gel electrophoresis of biotinylated cell lysates (PAGE) followed by mass spectrometry (MS) and antibody array analysis (MAP).

Liang R, Arif T, Kalmykova S et al. Restraining Lysosomal Activity Preserves Hematopoietic Stem Cell Quiescence and Potency Cell Stem Cell 2020-03-05 [PMID: 32109377]

Ahonen M, Haridas P, Mysore R et al. miR-107 inhibits CDK6 expression, differentiation, and lipid storage in human adipocytes. Mol. Cell. Endocrinol. 2018-09-24 [PMID: 30261211] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-87262

NBP1-87262PEP	Cdk6 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-87262

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

