

Product Datasheet

Indoleamine 2,3-dioxygenase/IDO Antibody - BSA Free NBP1-87703

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-87703

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-87703



NBP1-87703**Indoleamine 2,3-dioxygenase/IDO Antibody - BSA Free****Product Information**

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Description	Novus Biologicals Rabbit Indoleamine 2,3-dioxygenase/IDO Antibody - BSA Free (NBP1-87703) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-Indoleamine 2,3-dioxygenase/IDO Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3620
Gene Symbol	IDO1
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: MAHAMENSWTISKEYHIDEVGFALPNPQENLPDFYNDWMFIAKHLPDLIESG QLRERVEKLNMLSIDHLTDHKSQRLARLVLCITMAYVWGKGHGDVRKVL

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:1000 - 1:2500, Immunohistochemistry-Paraffin 1:1000 - 1:2500, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Western Blot: Indoleamine 2,3-dioxygenase/IDO Antibody [NBP1-87703]
- Analysis in EFO-21 cells transfected with control siRNA, target specific siRNA probe #1 and #2, using anti-IDO1 antibody. Remaining relative intensity is presented. Loading control: anti-PPIB.




Immunohistochemistry-Paraffin: Indoleamine 2,3-dioxygenase/IDO Antibody [NBP1-87703] - Staining in human placenta and cerebral cortex tissues . Corresponding IDO1 RNA-seq data are presented for the same tissues.




Western Blot: Indoleamine 2,3-dioxygenase/IDO Antibody [NBP1-87703] - Analysis in control (vector only transfected HEK293T lysate) and IDO1 over-expression lysate (Co-expressed with a C-terminal myc-DDK tag (3.1 kDa) in mammalian HEK293T cells).




Immunohistochemistry-Paraffin: Indoleamine 2,3-dioxygenase/IDO Antibody [NBP1-87703] - Staining of human placenta shows moderate to strong cytoplasmic positivity in endothelial cells.




Immunohistochemistry-Paraffin: Indoleamine 2,3-dioxygenase/IDO Antibody [NBP1-87703] - Staining of human cerebral cortex shows no cytoplasmic positivity in neuronal cells as expected.



Immunohistochemistry-Paraffin: Indoleamine 2,3-dioxygenase/IDO Antibody [NBP1-87703] - Staining of human lymph node shows strong positivity in non - germinal center cells.



Immunohistochemistry-Paraffin: Indoleamine 2,3-dioxygenase/IDO Antibody [NBP1-87703] - Staining of human small intestine shows moderate to strong positivity in lymphoid cells.



Publications

Workman MJ, Troisi E, Targan SR, et al. Modeling Intestinal Epithelial Response to Interferon- γ in Induced Pluripotent Stem Cell-Derived Human Intestinal Organoids International journal of molecular sciences 2020-12-30 [PMID: 33396621] (WB, ICC/IF, Human)

Mehra S, Alvarez X, Didier PJ et al. Granuloma Correlates of Protection Against Tuberculosis and Mechanisms of Immune Modulation by Mycobacterium tuberculosis. J Infect Dis 2013-04-01 [PMID: 23255564]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-87703

NBP1-87703PEP	Indoleamine 2,3-dioxygenase/IDO Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-87703

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

