

Product Datasheet

Kininogen Antibody - BSA Free NBP1-89993

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-89993

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-89993



NBP1-89993

Kininogen Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Rabbit
Gene ID	3827
Gene Symbol	KNG1
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: LFLTPDCKSLWNGDTGECTDNAYIDIQLRIASFQNCDIYPGKDFVQPPTKICVG CPRDIPTNSPELEETLTHITITKLNAENNATFYFKIDNVKKARVQVVAGKKYFIDFV ARETTCSKESNEELTESCETK

Product Application Details


Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:500 - 1:1000, Immunohistochemistry-Paraffin 1:500 - 1:1000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Western Blot: Kininogen Antibody [NBP1-89993] - Analysis in control (vector only transfected HEK293T lysate) and Kininogen over-expression lysate (Co-expressed with a C-terminal myc-DDK tag (3.1 kDa) in mammalian HEK293T cells).




Immunohistochemistry-Paraffin: Kininogen Antibody [NBP1-89993] - Staining of human kidney using Kininogen antibody NBP1-89993.




Western Blot: Kininogen Antibody [NBP1-89993] - Lane 1: Marker [kDa] 250, 130, 95, 72, 55, 36, 28, 17, 10. Lane 2: Negative control (vector only transfected HEK293T lysate). Lane 3: Over-expression lysate (Co-expressed with a C-terminal myc-DDK tag (3.1 kDa) in mammalian HEK293T cells).




Immunohistochemistry-Paraffin: Kininogen Antibody [NBP1-89993] - Staining of human skeletal muscle shows no positivity in myocytes as expected.




Immunohistochemistry-Paraffin: Kininogen Antibody [NBP1-89993] - Staining in human kidney and skeletal muscle tissues.. Corresponding Kininogen RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: Kininogen Antibody [NBP1-89993] - Staining of human colon, kidney, prostate and testis using Kininogen antibody NBP1-89993 (A) shows similar protein distribution across tissues to independent antibody NBP1-89994 (B).




Immunohistochemistry-Paraffin: Kininogen Antibody [NBP1-89993] - Staining of human prostate.



Immunohistochemistry-Paraffin: Kininogen Antibody [NBP1-89993] - Staining of human colon.



Immunohistochemistry-Paraffin: Kininogen Antibody [NBP1-89993] - Staining of human testis.



Publications

Kampf C, Mardinoglu A, Fagerberg L et al. The human liver-specific proteome defined by transcriptomics and antibody-based profiling. FASEB J 2014-07-01 [PMID: 24648543]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-89993

NBP1-89993PEP	Kininogen Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-89993

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

