

Product Datasheet

TAF15 Antibody (4D71) [PE]

NBP1-92685PE

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-92685PE

Updated 10/29/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-92685PE



NBP1-92685PE

TAF15 Antibody (4D71) [PE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	4D71
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	PE
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	68 kDa
Product Description	
Description	This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Host	Mouse
Gene ID	8148
Gene Symbol	TAF15
Species	Human, Mouse, Rat
Immunogen	Full length recombinant human TAF15 expressed in and purified from E. coli. [UniProt# Q92804]
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence
Application Notes	This TAF15 (4D71) antibody is useful for Immunocytochemistry/Immunofluorescence and Western blot, where a band can be seen at approximately 68 kDa. The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-92685PE

NBP1-97005PE	Mouse IgG1 Isotype Control (MG1) [PE]
NBP1-92685F	TAF15 Antibody (4D71) [FITC]
NBP2-22848	Recombinant Human TAF15 His Protein
H00006908-Q01-10ug	Recombinant Human TATA binding protein TBP GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-92685PE

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

