

# Product Datasheet

## Recombinant Human HOMER3 His Protein NBP1-98979-100ug

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-98979](http://www.novusbio.com/NBP1-98979)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-98979](http://www.novusbio.com/reviews/destination/NBP1-98979)



**NBP1-98979-100ug**

Recombinant Human HOMER3 His Protein

**Product Information**

|                                |                                                                                        |
|--------------------------------|----------------------------------------------------------------------------------------|
| <b>Unit Size</b>               | 100 ug                                                                                 |
| <b>Concentration</b>           | 0.25 mg/ml                                                                             |
| <b>Storage</b>                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| <b>Preservative</b>            | No Preservative                                                                        |
| <b>Purity</b>                  | >90%, by SDS-PAGE                                                                      |
| <b>Buffer</b>                  | 20 mM Tris-HCl buffer (pH8.0), 30% glycerol, 0.15 M NaCl, 1 mM DTT                     |
| <b>Target Molecular Weight</b> | 42.5 kDa                                                                               |

**Product Description**


|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | <p>A recombinant protein with a N-Terminal His-tag and corresponding to the amino acids 1-361 of Human HOMER3</p> <p><b>Source:</b> <i>E.coli</i></p> <p><b>Amino Acid Sequence:</b> MGSSHHHHHH SSGLVPRGSH MGSHMMSTAR EQPIFSTRAH VFQIDPATKR NWIPAGKHAL TVSYFYDATR NVYRIISIGG AKAIINSTVT PNMTFTKTSQ KFGQWADSRA NTVYGLGFAS EQHLTQFAEK FQEVKEAARL AREKSQDGGE LTSPALGLAS HQVPPSPLVS ANGPGEELKF RSQSADAPGP TERERLKKML SEGSGVEVQW EAEFFALQDS NNKLAGALRE ANAAAAQWRQ QLEAQRAEAE RLRQRVAELE AQAASEVTPT GEKEGLGQGQ SLEQLEALVQ TKDQEIQTLK SQTGGPREAL EAAEREETQQ KVQDLETRNA ELEHQLRAME RSLEEARAER ERARA EVGRA AQLLDVSLFE LSELREGLAR LAEAAP</p> |
| <b>Gene ID</b>     | 9454                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Gene Symbol</b> | HOMER3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Species</b>     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |

**Product Application Details**

|                              |          |
|------------------------------|----------|
| <b>Applications</b>          | SDS-Page |
| <b>Recommended Dilutions</b> | SDS-Page |

**Images**

SDS-Page: HOMER3 Protein [NBP1-98979]





### Novus Biologicals USA

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### Bio-Techne Canada

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### Bio-Techne Ltd

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### General Contact Information

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### Products Related to NBP1-98979-100ug

|               |                                           |
|---------------|-------------------------------------------|
| NBP2-32607PEP | HOMER3 Recombinant Protein Antigen        |
| NBP2-13075    | beta Amyloid Antibody (MOAB-2) - BSA Free |
| NBP2-32607    | HOMER3 Antibody                           |
| NB100-1965    | Calnexin Antibody - BSA Free              |

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-98979](http://www.novusbio.com/reviews/submit/NBP1-98979)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

