

Product Datasheet

UNG Antibody (OTI2C12) NBP2-01230

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-01230

Updated 9/6/2018 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-01230



NBP2-01230

UNG Antibody (OTI2C12)

Product Information

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI2C12
Preservative	0.02% Sodium Azide
Isotype	IgG2b
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	34.5 kDa

Product Description


Host	Mouse
Gene ID	7374
Gene Symbol	UNG
Species	Human
Reactivity Notes	Human
Immunogen	Full length recombinant protein of human UNG (NP_550433) produced in HEK293T cells.

Product Application Details


Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-2000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:150

Images


Western Blot: uracil-DNA glycosylase Antibody (2C12) [NBP2-01230]
 Analysis of extracts (35ug) from 9 different cell lines by using anti-uracil-DNA glycosylase monoclonal antibody.




Immunocytochemistry/Immunofluorescence: uracil-DNA glycosylase Antibody (2C12) [NBP2-01230] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY uracil-DNA glycosylase.




Immunohistochemistry-Paraffin: uracil-DNA glycosylase Antibody (2C12) [NBP2-01230] - Staining of paraffin-embedded Human liver tissue using anti-uracil-DNA glycosylase mouse monoclonal antibody.



Western Blot: uracil-DNA glycosylase Antibody (2C12) [NBP2-01230] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY uracil-DNA glycosylase (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-uracil-DNA glycosylase.



Immunohistochemistry-Paraffin: uracil-DNA glycosylase Antibody (2C12) [NBP2-01230] - Staining of paraffin-embedded Human Kidney tissue using anti-uracil-DNA glycosylase mouse monoclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-01230

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

