

Product Datasheet

RALBP1 Antibody (OTI11B2)

NBP2-02293

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-02293

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-02293



NBP2-02293**RALBP1 Antibody (OTI11B2)**


Product Information	
Unit Size	0.1 ml
Concentration	0.84 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI11B2
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	76 kDa

Product Description	
Host	Mouse
Gene ID	10928
Gene Symbol	RALBP1
Species	Human, Mouse, Rat, Canine
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Immunogen	Full length human recombinant protein of human RALBP1 (NP_006779) produced in HEK293T cell.


Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2000, Flow Cytometry 1:100, Immunohistochemistry 1:50, Immunocytochemistry/ Immunofluorescence 1:100, Immunoprecipitation 4ug/ml, Immunohistochemistry-Paraffin 1:50

Images


Western Blot: RALBP1 Antibody (OTI11B2) [NBP2-02293] - Analysis of extracts (35ug) from 9 different cell lines by using anti-RALBP1 monoclonal antibody.




Immunocytochemistry/Immunofluorescence: RALBP1 Antibody (OTI11B2) [NBP2-02293] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY RALBP1.




Immunohistochemistry-Paraffin: RALBP1 Antibody (OTI11B2) [NBP2-02293] - Staining of paraffin-embedded prostate tissue using anti-RALBP1 mouse monoclonal antibody.




Flow Cytometry: RALBP1 Antibody (OTI11B2) [NBP2-02293] - Analysis of Jurkat cells, using anti-RALBP1 antibody, (Red) compared to a nonspecific negative control antibody (Blue).




Western Blot: RALBP1 Antibody (OTI11B2) [NBP2-02293] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RALBP1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RALBP1.




Immunohistochemistry-Paraffin: RALBP1 Antibody (OTI11B2) [NBP2-02293] - Staining of paraffin-embedded breast tissue using anti-RALBP1 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: RALBP1 Antibody (OTI11B2) [NBP2-02293] - Staining of paraffin-embedded Carcinoma of bladder tissue using anti-RALBP1 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: RALBP1 Antibody (OTI11B2) [NBP2-02293] - Staining of paraffin-embedded Carcinoma of prostate tissue using anti-RALBP1 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: RALBP1 Antibody (OTI11B2) [NBP2-02293] - Staining of paraffin-embedded endometrium tissue using anti-RALBP1 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: RALBP1 Antibody (OTI11B2) [NBP2-02293] - Staining of paraffin-embedded Ovary tissue using anti-RALBP1 mouse monoclonal antibody.



Immunoprecipitation: RALBP1 Antibody (OTI11B2) [NBP2-02293] - (Negative control: IP without adding anti-RALBP1 antibody.). For each experiment, 500ul of DDK tagged RALBP1 overexpression lysates (at 1:5 dilution with HEK293T lysate), 2ug of anti-RALBP1 antibody and 20ul (0.1mg) of goat anti-mouse conjugated magnetic beads were mixed and incubated overnight. After extensive wash to remove any non-specific binding, the immuno-precipitated products were analyzed with rabbit anti-DDK polyclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-02293

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP2-30698PEP	RALBP1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-02293

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

