

Product Datasheet

DNAJA2 Antibody (OTI3A10)

NBP2-02996

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-02996

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-02996



NBP2-02996

DNAJA2 Antibody (OTI3A10)

Product Information

Unit Size	0.1 ml
Concentration	0.37 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI3A10
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	45.6 kDa

Product Description


Host	Mouse
Gene ID	10294
Gene Symbol	DNAJA2
Species	Human, Mouse, Rat, Canine, Monkey
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Immunogen	Full length human recombinant protein of human DNAJA2 (NP_005871) produced in HEK293T cell.

Product Application Details

Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500, Flow Cytometry 1:100, Immunohistochemistry 1:50-150, Immunocytochemistry/ Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:50-150

Images


Western Blot: DNAJA2 Antibody (3A10) [NBP2-02996] Analysis of extracts (35ug) from 9 different cell lines by using anti-DNAJA2 monoclonal antibody.




Immunocytochemistry/Immunofluorescence: DNAJA2 Antibody (3A10) [NBP2-02996] - Staining of COS7 cells transiently transfected by pCMV6-ENTRY DNAJA2.




Immunohistochemistry-Paraffin: DNAJA2 Antibody (3A10) [NBP2-02996] - Staining of paraffin-embedded Human thyroid tissue using anti-DNAJA2 mouse monoclonal antibody.




Flow Cytometry: DNAJA2 Antibody (3A10) [NBP2-02996] - HEK293T cells transfected with either overexpression plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-DNAJA2 antibody, and then analyzed by flow cytometry.




Western Blot: DNAJA2 Antibody (3A10) [NBP2-02996] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DNAJA2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DNAJA2.




Immunohistochemistry-Paraffin: DNAJA2 Antibody (3A10) [NBP2-02996]
- Staining of paraffin-embedded Adenocarcinoma of Human breast tissue using anti-DNAJA2 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: DNAJA2 Antibody (3A10) [NBP2-02996]
- Staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-DNAJA2 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: DNAJA2 Antibody (3A10) [NBP2-02996]
- Staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-DNAJA2 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: DNAJA2 Antibody (3A10) [NBP2-02996]
- Staining of paraffin-embedded Carcinoma of Human thyroid tissue using anti-DNAJA2 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: DNAJA2 Antibody (3A10) [NBP2-02996]
- Staining of paraffin-embedded Human breast tissue using anti-DNAJA2 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: DNAJA2 Antibody (3A10) [NBP2-02996]
- Staining of paraffin-embedded Human Kidney tissue using anti-DNAJA2 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: DNAJA2 Antibody (3A10) [NBP2-02996]
- Staining of paraffin-embedded Human lymph node tissue using anti-DNAJA2 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: DNAJA2 Antibody (3A10) [NBP2-02996]
- Staining of paraffin-embedded Human lymphoma tissue using anti-DNAJA2 mouse monoclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-02996

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
H00010294-P01-10ug	Recombinant Human DNAJA2 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-02996

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

