

# Product Datasheet

## **GITR/TNFRSF18 Antibody (RM0066-7C13) - Azide and BSA Free** **NBP2-12338**

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-12338](http://www.novusbio.com/NBP2-12338)

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-12338](http://www.novusbio.com/reviews/destination/NBP2-12338)



**NBP2-12338**

GITR/TNFRSF18 Antibody (RM0066-7C13) - Azide and BSA Free

**Product Information**

<b>Unit Size</b>	0.1 mg
<b>Concentration</b>	LYOPH mg/ml
<b>Storage</b>	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	RM0066-7C13
<b>Preservative</b>	No Preservative
<b>Reconstitution Instructions</b>	Reconstitute with sterilized PBS to a final concentration of 0.5 mg/ml.
<b>Isotype</b>	IgG2
<b>Purity</b>	Protein A or G purified
<b>Buffer</b>	Lyophilized from a 0.2 um filtered solution in PBS. 0.025 mg size is provided in liquid form, PBS

**Product Description**


<b>Host</b>	Rat
<b>Gene ID</b>	8784
<b>Gene Symbol</b>	TNFRSF18
<b>Species</b>	Mouse
<b>Immunogen</b>	Mouse T cells

**Product Application Details**

<b>Applications</b>	Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready
<b>Recommended Dilutions</b>	Flow Cytometry 1:50-1000, Immunohistochemistry 1:50-200, Immunohistochemistry-Paraffin 1:50-1:200, CyTOF-ready
<b>Application Notes</b>	This antibody is Cytot ready.

**Images**

Immunohistochemistry-Paraffin: GITR/TNFRSF18 Antibody (RM0066-7C13) [NBP2-12338] - 4% PFA fixed/Paraffin embedded mouse Spleen, Lung and Kidney sections from LPS treated mice was subjected to immunohistochemistry staining (DAB) of GITR.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-12338**

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
NBP1-75398	Goat anti-Rat IgG (H+L) Secondary Antibody (Pre-adsorbed)
NBP2-12338UV	GITR/TNFRSF18 Antibody (RM0066-7C13) [DyLight 350]
NBP2-52215-0.05mg	Recombinant Human GITR/TNFRSF18 His Protein

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-12338](http://www.novusbio.com/reviews/submit/NBP2-12338)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

