

Product Datasheet

BAG3 Antibody NBP2-15543

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-15543

Updated 10/23/2024 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-15543




NBP2-15543**BAG3 Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 10% Glycerol
Target Molecular Weight	62 kDa
Product Description	
Host	Rabbit
Gene ID	9531
Gene Symbol	BAG3
Species	Human, Mouse
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human BAG3. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunoblotting, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Immunoblotting 1:100 - 1:2000




Images


Western Blot: BAG3 Antibody [NBP2-15543] - BAG3 antibody detects BAG3 protein by Western blot analysis. A. 50 ug mouse heart lysate/extract 7. 5 % SDS-PAGEBAG3 antibody dilution: 1:1000.




Immunohistochemistry-Paraffin: BAG3 Antibody [NBP2-15543] - BAG3 in high-grade serous carcinoma. 1:100 dilution, pH=6 retrieval buffer, Dako FLEX+ protocol. This image was submitted via customer Review.



Immunohistochemistry-Paraffin: BAG3 Antibody [NBP2-15543] - Immunohistochemical analysis of paraffin-embedded Cal27 xenograft, using antibody at 1:100 dilution.



Western Blot: BAG3 Antibody [NBP2-15543] - Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with BAG3 antibody (NBP2-15543) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Publications

Nymoén DA, Hetland Falkenthal TE, Holth A et al. Expression and clinical role of chemoresponse-associated genes in ovarian serous carcinoma. Gynecol. Oncol. 2015-07-29 [PMID: 26232338] (IHC-P, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-15543

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-72276-50ug	Recombinant Human BAG3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-15543

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

