

Product Datasheet

ARC/NOL3 Antibody NBP2-19561

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-19561

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-19561



NBP2-19561

ARC/NOL3 Antibody

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol, 1% BSA
Target Molecular Weight	24 kDa

Product Description


Host	Rabbit
Gene ID	8996
Gene Symbol	NOL3
Species	Human, Mouse, Rat
Immunogen	Recombinant protein encompassing a sequence within the center region of human ARC/NOL3. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Immunohistochemistry-Frozen 1:100-1:1000

Images


Immunocytochemistry/Immunofluorescence: ARC/NOL3 Antibody [NBP2-19561] - Confocal immunofluorescence analysis (Olympus FV10i) of paraformaldehyde-fixed HeLa, using NOL3 antibody (Green) at 1:500 dilution. Alpha-tubulin filaments were labeled with (Red) at 1:2000.




Immunohistochemistry-Paraffin: ARC/NOL3 Antibody [NBP2-19561] - Mouse brain. NOL3 stained by NOL3 antibody diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min.




Immunohistochemistry-Paraffin: ARC/NOL3 Antibody [NBP2-19561] - Hela xenograft, using NOL3 antibody at 1:500 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Immunohistochemistry-Paraffin: ARC/NOL3 Antibody [NBP2-19561] - NOL3 antibody detects NOL3 protein at cytoplasm and nucleus by immunohistochemical analysis. Sample: Paraffin-embedded rat brain. NOL3 stained by NOL3 antibody diluted at 1:500.




Immunohistochemistry-Frozen: ARC/NOL3 Antibody [NBP2-19561] - NOL3 antibody detects NOL3 protein by immunohistochemical analysis. Sample: Frozen-sectioned mouse cerebellum. Green: NOL3 stained by NOL3 antibody diluted at 1:250. Red: NF-H, stained by NF-H antibody diluted at 1:500. Blue: Fluoroshield with DAPI.




Western Blot: ARC/NOL3 Antibody [NBP2-19561] - Various whole cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with ARC/NOL3 antibody (NBP2-19561) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Immunohistochemistry-Paraffin: ARC/NOL3 Antibody [NBP2-19561] - ARC/NOL3 antibody detects ARC/NOL3 protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded mouse brain. ARC/NOL3 stained by ARC/NOL3 antibody (NBP2-19561) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: ARC/NOL3 Antibody [NBP2-19561] - ARC/NOL3 antibody detects ARC/NOL3 protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded rat brain. ARC/NOL3 stained by ARC/NOL3 antibody (NBP2-19561) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-19561

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-23300	Recombinant Human ARC/NOL3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-19561

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

