

Product Datasheet

Nucleoporin 107 Antibody - BSA Free NBP2-19602

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-19602

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-19602



NBP2-19602**Nucleoporin 107 Antibody - BSA Free**

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 20% Glycerol
Target Molecular Weight	106 kDa


Product Description	
Description	Novus Biologicals Rabbit Nucleoporin 107 Antibody - BSA Free (NBP2-19602) is a polyclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	57122
Gene Symbol	NUP107
Species	Human
Reactivity Notes	Zebrafish (80%), Xenopus laevis (84%), Chicken (87%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human Nucleoporin 107. The exact sequence is proprietary.

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:3000, Immunoprecipitation 1:100-1:500




Images

Western Blot: Nucleoporin 107 Antibody [NBP2-19602] - Sample (30 ug of whole cell lysate) A: Hela 7. 5% SDS PAGE gel, diluted at 1:1000.



Immunoprecipitation: Nucleoporin 107 Antibody [NBP2-19602] - A : 30 ug whole cell lysate/extract of NUP107 protein expressing 293T cells B : Control with 3 ug of pre-immune rabbit IgG C : Immunoprecipitation of NUP107 by 3 ug of NUP107 antibody . 7.5% SDS-PAGE The immunoprecipitated NUP107 protein was detected by NUP107 antibody diluted at 1 : 1000.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-19602

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-76926PEP	Nucleoporin 107 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-19602

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

