

Product Datasheet

PAK6 Antibody NBP2-19685

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-19685

Updated 1/5/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-19685




NBP2-19685**PAK6 Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol
Target Molecular Weight	75 kDa
Product Description	
Host	Rabbit
Gene ID	56924
Gene Symbol	PAK6
Species	Human
Immunogen	Recombinant protein encompassing a sequence within the center region of human PAK6. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000-1:10000, Immunohistochemistry 1:100-1:1000, Immunoprecipitation 1:100-1:500, Immunohistochemistry-Paraffin 1:100-1:1000




Images


Western Blot: PAK6 Antibody [NBP2-19685] - Sample (30 ug of whole cell lysate) A: Non-transfected 293T lysates B: PAK6 transfected 293T lysates 7. 5% SDS PAGE gel, diluted at 1:5000.



Immunohistochemistry-Paraffin: PAK6 Antibody [NBP2-19685] - Immunohistochemical analysis of paraffin-embedded HeLa xenograft, using antibody at 1:500 dilution.



Immunoprecipitation: PAK6 Antibody [NBP2-19685] - PC-3 whole cell extract A. Control with 4 ug of preimmune Rabbit IgG B. Immunoprecipitation of PAK6 protein by 4 ug PAK6 antibody 7.5 % SDS-PAGE The immunoprecipitated PAK6 protein was detected by PAK6 antibody diluted at 1:500. [EasyBlot anti-rabbit IgG was used as a secondary reagent]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-19685

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

