

Product Datasheet

PGLS Antibody NBP2-19785

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-19785

Updated 3/16/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-19785




NBP2-19785**PGLS Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1% BSA, 20% Glycerol
Target Molecular Weight	28 kDa
Product Description	
Host	Rabbit
Gene ID	25796
Gene Symbol	PGLS
Species	Human, Mouse
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Rat (88%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human PGLS. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:1000-1:10000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000




Images


Western Blot: PGLS Antibody [NBP2-19785] - Sample (30 ug of whole cell lysate) A: SK-N-SH 12% SDS PAGE gel, diluted at 1:5000.




Immunocytochemistry/Immunofluorescence: PGLS Antibody [NBP2-19785] - HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: PGLS protein stained by PGLS antibody diluted at 1:500. Blue: Hoechst 33342 staining. Scale bar = 10 μ m.



Immunohistochemistry-Paraffin: PGLS Antibody [NBP2-19785] - : Paraffin-embedded mouse prostate. PGLS antibody dilution: 1:500.



Immunohistochemistry-Paraffin: PGLS Antibody [NBP2-19785] - Immunohistochemical analysis of paraffin-embedded U87 xenograft, using antibody at 1:500 dilution.



Publications

Pushpak Bhattacharjee, Die Wang, Dovile Anderson, Joshua N Buckler, Eveline de Geus, Feng Yan, Galina Polekhina, Ralf Schittenhelm, Darren J Creek, Lawrence D Harris, Anthony J Sadler The immune response to RNA suppresses nucleic acid synthesis by limiting ribose 5-phosphate The EMBO Journal 2024-05-22 [PMID: 38778156]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-19785

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-19785

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

